

Analog Adapter

PS-2158



Introduction

The PS-2158 Analog Adapter allows an analog ScienceWorkshop sensor to be used with a PASPORT interface. Multiple Analog Adapters can be used simultaneously to support multiple ScienceWorkshop sensors.

| Compatible PASPORT Interfaces | Part Number |
|-------------------------------|-------------|
| Xplorer GLX ¹ | PS-2002 |
| Xplorer ^{1, 2} | PS-2000 |
| Power Link | PS-2001 |
| USB Link | PS-2100 |

¹When the Analog Adapter is used with an Xplorer or Xplorer GLX in logging mode (not connected to a computer), the logger records output voltage from the sensor. A future free upgrade to the GLX will support calibrated data logging. Visit www.pasco.com for the latest GLX upgrade.

²When the Analog Adapter is used with an Xplorer (PS-2000) in logging mode, the Xplorer Power Adapter (PS-2530) is recommended (and in some cases required) due to the high current draw of some ScienceWorkshop sensors.

Equipment and Software Setup

DataStudio version 1.9.5 (or later) is required. Visit www.pasco.com to download the latest version of DataStudio.

The first three steps can be performed in any order.

- **1.** Connect an analog ScienceWorkshop sensor to the cable of the Analog Adapter.
- 2. Connect the Analog Adapter to a PASPORT interface.
- **3.** Connected the PASPORT interface to a computer.
- 4. If DataStudio is not running, it will launch automatically.
- **5.** The ScienceWorkshop Analog Sensors list will open automatically. Select the connected sensor from the list and click OK.
- 6. To view and change the sampling rate, gain, and other settings of the sensor and adapter, press the Setup button in DataStudio's main tool bar; the Experiment Setup window will open. Refer to DataStudio documentation for further instructions.

|:/_{_\}+\) • °

800-772-8700 www.pasco.com

Specifications

| Maximum input voltage range ¹ | -10 V to +10 V |
|---|--|
| Absolute maximum input voltage range without damage | -40 V to +40 V |
| Input impedance | 1 MΩ |
| Gain | 1, 10, and 100 (selectable in DataStudio) |
| Maximum Sampling Rate | 50 kHz with Xplorer GLX 1000 Hz with other interfaces |
| Analog-to-Digital Conversion | 12 bit |
| Voltage Resolution | 5 mV at gain = 1 0.5 mV at gain = 10 0.05 mV at gain = 100 |
| Offset voltage accuracy | < ±5 mV |
| Full-scale voltage accuracy | < ±15 mV |

¹When used with a ScienceWorkshop Voltage Sensor (CI-6503) or Current Sensor (CI-6556), the sensor measures differential voltage, but the maximum input voltage range to either terminal is ± 10 V.

Technical Support

For assistance with any PASCO product, contact PASCO at:

| Address: | PASCO scientific |
|----------|--------------------------|
| | 10101 Foothills Blvd. |
| | Roseville, CA 95747-7100 |
| Phone: | 916-786-3800 (worldwide) |
| | 800-772-8700 (U.S.) |
| Fax: | (916) 786-3292 |
| Web: | www.pasco.com |
| Email: | support@pasco.com |

Limited Warranty

For a description of the product warranty, see the PASCO catalog.

Copyright

The PASCO scientific 012-08921A Analog Adapter Instruction Sheet is copyrighted with all rights reserved. Permission is granted to non-profit educational institutions for reproduction of any part of this manual, providing the reproductions are used only in their laboratories and classrooms, and are not sold for profit. Reproduction under any other circumstances, without the written consent of PASCO scientific, is prohibited.

Trademarks

PASCO, PASCO scientific, DataStudio, PASPORT, and ScienceWorkshop are trademarks or registered trademarks of PASCO scientific, in the United States and/or in other countries. All other brands, products, or service names are or may be trademarks or service marks of, and are used to identify, products or services of, their respective owners. For more information visit www.pasco.com/legal.