**USB Interface**  
Instruction Manual

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## INTRODUCTION

This Interface Kit is for use with the following Ohaus products: Scout Pro, Traveler and Navigator.

The Ohaus USB Interface Kit is a unique solution to the problem of connecting a balance to a computer using a Universal Serial Bus (USB). USB devices are categorized into classes such as disk drives, digital cameras, printers, etc. Balances do not have a commonly used class so the Ohaus USB interface uses a generic interface based on the RS232 serial standard.

Data sent from the balance to a computer is in USB format. The USB data is directed to a *virtual port*. This port then appears as an RS232 port to the application program.

When sending a command from a computer to the balance, the application program sends a command to the *virtual port* as if it were an RS232 port. The computer then directs the command from the *virtual port* to the computers USB connector where the balance is connected. The Adapter Kit receives the USB signal and reacts to the command.

The USB Interface kit includes a CD with the software drivers to create the required *virtual port* on the computer.


## SYSTEM REQUIREMENTS

- PC running Windows 98®, Windows 98SE®, Windows ME®, Windows 2000®, Windows XP® or Windows 7®
- Available USB port (Type A, 4-pin, female)


## INTERFACE INSTALLATION

Install the Interface module on the underside of the balance as shown.


Note: The appearance of your model may be different.


1. Slide Blank Cover to rear.

2. Remove the Blank Cover.

3. Install Interface Module.

4. Slide Module forward.

## USB CONNECTION

The USB cable terminates with a 4-pin, male, USB Type A connector.

1. Power on the computer and verify that the USB port is enabled and working properly.
2. Plug the cable's USB connector into the computer's USB connector. Windows® should detect a USB device and the New Hardware Wizard will be initialized.

## VIRTUAL PORT SOFTWARE INSTALLATION

1. Insert the supplied CD into the CD drive.

Different versions of Windows® have slightly different steps to load the driver that is on the CD. In all versions the New Hardware Wizard guides you through the required steps to select the driver that is located on the CD.

Example of Windows XP Hardware Wizard

2. After clicking Finish, the virtual port should be ready for use. Windows® typically adds the virtual port in sequence after the highest number COM port. For example, on PC's equipped with up to 4 COM ports, the virtual port will be COM5. When using the USB interface with programs that limit the number of COM port designations (e.g. Ohaus MassTracker allows only COM1, 2, 3, & 4), it may be necessary to assign one of these port numbers to the new virtual port. This can be done in the Port Settings of the Device Manager utility, found in the Windows Control Panel.

## BALANCE SETUP

- A) The Ohaus Interface is preset to communicate using the following settings: 2400 baud, 7 bit, no parity, no handshake. If different settings are desired, it will be necessary to change either the balance settings, or the computer settings.
- B) Upon installation of the interface, the balance will recognize the USB Interface and add relevant items to the menu. Configure the balance to the desired USB and printing parameters; refer to the balance User Manual for assistance in using the menus.

### PRINT

Stable  
 On, Off  
 A-Print  
 Cont, On.Stbl, On.Acc\*, 5sec, 15sec,  
 30sec, 60sec, off  
 End  
 Yes, No

\*Note: Print and USB menu selections may vary between the different balance series.

### USB

USB  
 On / off  
 Baud  
 600, 1200, 2400, 4800, 9600, 19200  
 Parity  
 7-even, 7-odd, 7-none, 8-none  
 Handshake  
 None, Xon-Xoff, RTS-CTS  
 END  
 Yes, No

### SETUP DEFINITIONS

PRINT / Stable - On  
 PRINT / Auto Print – Continuous  
 PRINT / Auto Print – On Stable  
 PRINT / Auto Print – (xx) sec  
 PRINT / Auto Print – Off  
 USB / USB –Off  
 USB / Baud, Parity, Handshake

Balance transmits stable data only.  
 Balance repeatedly sends data as fast as possible.  
 Balance automatically outputs data upon stability.  
 Balance sends data every (xx) seconds.  
 Only sends data when PRINT is pressed.  
 Turns Interface off, saving battery power.  
 Set parameters to match printer or computer.

## OPERATION

### USB OUTPUT

With the Adapter installed the balance will operate in one of the following three ways:

- USB = on, Auto Print = off, Stable = on or off  
 Pressing PRINT will send the display data based on the settings in the Menu.  
 If Stable is set to On, the balance will wait for a stable reading before sending data.
- USB = on, Auto Print = on, Stable = on or off  
 The balance will automatically send data.  
 If Stable is set to On, only stable values will be sent.
- USB = off  
 The interface is turned off. During battery operation, this saves significant power.

Data sent from the Interface is in standard ASCII format terminated with a [CRLF] (carriage return-line feed). The output format is as follows:

**Scout Pro (SPx) Output:****Output Format 1 (for the models 303/123/202/402/602/2001/6001/401FZH/601FZH/6000FZH):**Output String:

[weight]	12 characters (right justified)
[space]	1 character
[unit]	5 characters
[space]	1 character
[stability indicator]	1 character; "?" when unstable, blank when stable
[space]	1 character
[Legend]	10 character

**Note:** All of the fields have fixed length.

Output examples (space indicated by '\*'):

```

*****0.00*g*****
*****12.73*g*****?
*****12.72*g*****?
*****0.4500*oz*****
*****0.5950*oz*****?
*****3*PCS**?
*****3*PCS***
*****0.85*oz*****WET*WT
*****0.85*oz*****AIR*WT

```

**Output Format 2 (for the models 401/601/6000):**Output String:

[weight]	11~12 characters (right justified)
[space]	1 character
[unit]	1~5 characters
[space]	1 character
[stability indicator]	1 character; "?" when unstable, blank when stable
[space]	1 character
[Legend]	10 character

**Note:** The unit field length varies with different units. The weight field could be 11 or 12, depending on if the weight string has an additional decimal point.

Output example (space indicated by '\*'):

```

*****200*g**
*****311*g*?
*****0.211*kg*?

```

**Navigator (NV/NVL/NVT) Output:**

Output String:

[weight]	10 characters (right justified)
[space]	1 character
[unit]	1~5 characters (left justified)
[space]	1 character
[stability indicator]	1 character; "?" when unstable, blank when stable
[space]	1 character
[NET]	0 or 3 character
[Legend]	0~16 character

**Note:** The unit field length varies with different units, the Legend field is only printed in certain modes, and the time string only printed at internal print.

Output example (space indicated by '\*'):

```

*****200*g***
*****15.06*g***
*****15*g***NET
*****124*g?*NET
*****15*g***NET**00:00:02
*****122*g?*NET**00:00:03
*****15*g***NET*UNDER
***5:10.75*lb:oz**?*NET*ACCEPT*00:00:05
***0:*0.50*lb:oz**?*NET*UNDER*00:00:06
    
```

**Traveler (TAx/SE) Output:**

Output String:

[weight]	11 characters (right justified)
[space]	1 character
[unit]	1~5 characters
[space]	1 character
[stability indicator]	1 character; "?" when unstable, blank when stable
[space]	1 character
[Legend]	10 character

**Note:** The unit field length varies with different units.

Output example (space indicated by '\*'):

```

*****-0.01*g*?
*****-0.04*g*?
*****4.20*g*?
*****34.91*g*?
    
```

USB INPUT

The balance will respond to various commands sent via the interface adapter. Terminate the following commands when with a [CR] or [CRLF].

Scout Pro and Traveler Commands

- ? print current mode
- 0A** turn Auto-print off
- SA** Auto-print on, prints on stability
- CA** Continuous Auto-print
- (n)A** Auto-print on 1 to 3600 second intervals (n = 1 to 3600)
- C** perform span calibration
- L** perform linearity calibration
- 0M** gram mode
- 1M** ounce mode
- 2M** troy ounce mode
- 3M** pennyweight mode
- 4M** parts counting mode
- 5M** pound mode
- T** tare balance, same as pressing **ON-ZERO**
- V** print software version
- (Esc)R** reset balance to factory defaults
- P** same as pressing **PRINT**
- LE** print last error code, i.e. [**Err 0**]
- 0S** print unstable data
- 1S** print stable data only

Navigator Commands

- P** same as pressing **Print**
- SP** print stable weight only
- IP** immediate print of displayed weight (stable or unstable)
- CP** Continuous print of weights
- SLP** Auto-print stable non-zero weight only
- SLZP** Auto-print stable non-zero weight and zero reading
- xP** Auto-print on 1 to 3600 second intervals (x = 1 to 3600)
- 0P** turns auto-print off
- PM** print current mode
- M** advance to the next enabled mode
- PU** print current unit
- U** advance to the next enabled unit
- T** same as pressing **Tare**
- Z** same as pressing **Zero**
- PV** print software version

AUTO-PRINT OPERATION

Once Auto-Print is activated in the menu, the balance will send data as required. To temporarily stop Auto-Printing, press the PRINT key. If there is data in the print buffer the printer will finish printing this data. A second press will resume Auto-Printing.

## ACCESSORIES

For a complete listing of Ohaus printers and other accessories, contact Ohaus Corporation or visit [www.Ohaus.com](http://www.Ohaus.com)

## COMPLIANCE

This accessory has been tested and complies with the approvals listed in the applicable Instruction Manual.

### Disposal

In conformance with the European Directive 2002/96 EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

For disposal instructions in Europe, refer to [www.ohaus.com/weee](http://www.ohaus.com/weee).

Thank you for your contribution to environmental protection.

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P/N 83032208 B © 2016 Ohaus Corporation, all rights reserved / todos los derechos reservados / tous droits réservés / Alle Rechte vorbehalten / tutti i diritti riservati / все права защищены.