MULTI-FUNCTION INCUBATOR Series

ENVIRO-GENIE™

INCUBATOR GENIE™

Instruction Manual

EN

ES

FR

DE

IT

Change History

#	Date	Version	Descriptions
1	2026.03.08	A	New release

Table of Contents

Introduction	2
Applicable Products	2
Safety Information	2
Package Contents	3
Intended Use	4
Installation	4
Specifications	4
Control Panel	
<i>ENVIRO-GENIE</i>	5
<i>INCUBATOR GENIE</i>	7
Mixing Modes Configuration	8
Operation Instruction	
<i>ENVIRO-GENIE</i>	11
<i>INCUBATOR GENIE</i>	15
Maintenance	16
Technical Support Information	16
Technical Data	16
Compliance	17

Introduction

This manual contains installation, operation, and maintenance instructions for OHAUS MULTI-FUNCTION INCUBATORS. Please read the manual completely before using.

Applicable Products

This instruction manual is applicable to the following products:

Product	Model Name
ENVIRO-GENIE	e-MI3C (SI-xxxx*)
INCUBATOR GENIE	e-MI3 (SI-xxxx*)

* SI-xxxx represents the original SI model name.

Safety Information

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

WARNING For a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided.

CAUTION For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.

ATTENTION For important information about the product. May lead to equipment damage if not avoided.

NOTE For useful information about the product.

Warning Symbols

General hazard

Electric shock hazard

Safety Precautions

CAUTION: Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the rotary platform or the incubator lid.

WARNING! Unit Exceeds 18kg (40lbs.) Use at least two people to lift or move this unit to avoid injury.

Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly.

DO NOT immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

WARNING! The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer.

WARNING! DO NOT use the unit in a hazardous atmosphere or with hazardous materials for which the unit was not designed.

CAUTION: Replace fuse with the correct type and value only.

CAUTION: Ensure the power is off before shifting operating modes.

CAUTION: To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplugging from the wall outlet.

Earth Ground – Protective Conductor Terminal

Alternating Current

Package Contents

- Main Product
- 3x Clip Plates
- 4x Bag Mounting Strips
- 1x Non-Skid Tray
- 1x Shelf

We recommend you retain the original packaging for 90 days in case you need to return the product for any reason to your distributor or Ohaus.

Intended Use

This instrument is intended for use in laboratories, pharmacies, schools, businesses and light industry. It must only be used for processing materials as described in these operating instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from OHAUS, is considered as not intended. This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the instrument is not used according to these operating instructions, the intended protection provided by the instrument may be impaired.

Installation

Upon receiving the unit, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred during transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, place the unit on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand the typical heat produced by the unit. Do not position the equipment such that it is difficult to disconnect the power cord during use. Always place the unit on a clean and sturdy work surface.

The unit is supplied with a power cord that should be plugged into a properly grounded outlet. If the supplied power cord does not meet your needs, please use an approved power cord that suits local codes and electric supply. Replacement of the cord must be made by a qualified electrician.


Specifications

Product	ENVIRO-GENIE	INCUBATOR GENIE
Control	Digital	
Display	LCD	LED Display
Accessories Included	Clip Plate for 10-13 mm tubes, Clip Plate for 15-17 mm tubes, Clip Plate for 28-30 mm tubes, 4x Bag Mounting Strips, Non-Skid Tray, Shelf	
Capacity	4.5 kg / 10 lb. centered on axis for rotating / rocking, 2x 1L for stirring	4.5 kg / 10 lb. centered on axis for rotating / rocking
Motion	Rotate, Rock, Stir, Stir Reverse, Shake with Accessories	Rotate, Rock, Shake with Accessories
Temperature Range	4°C - 75°C	28°C - 75°C
Set Temperature Increments	0.1°C	1.0°C
Temperature Stability ($\pm$)	0.1°C	
Temperature accuracy ($\pm$)	0.2°C	
Uniformity	0.5°C	
Rotating Speed	1 rpm - 35 rpm	3 rpm - 35 rpm
Rocking Speed	2 rpm - 70 rpm	6 rpm - 70 rpm

Product	ENVIRO-GENIE	INCUBATOR GENIE
Rotating / Rocking Speed Increment	1 rpm	
Stirring Speed	1 rpm - 1000rpm	--
Stirring Speed Increment	1 rpm for speed 1-10 rpm 10 rpm for speed 10-1000 rpm	--
Timer	1 minute - 99 hours, 59 minutes	--
Timer Increment	1 minute	--
Stirring Reverse Timer	1 second – 60 seconds	--
Stirring Reverse Timer Increment	1 second	--
Alarm	Visible & audible, temperature & time	Visible & audible, temperature
Calibration	Temperature calibration	
Communication	RS232	--
Chamber Volume	24 L / 0.86 ft ³	
Inner Dimensions (L x W x H)	254 x 362 x 260 mm / 10 x 14.25 x 10.25 inch	
Tray Size (L x W)	203 x 305 mm / 8 x 12 inch	
Power	120V, 5A, 50/60Hz 230V, 2.5A, 50/60Hz 100V, 5A, 50/60Hz	
Dimensions (L x W x H)	495 x 622 x 368 mm / 19.5 x 24.5 x 14.5 inch	400 x 560 x 368 mm / 15.75 x 22 x 14.5 inch
Net Weight	36 kg / 80 lb.	19.5 kg / 43 lb.

Control Panel

ENVIRO-GENIE

Button Function

UP and **DOWN** buttons are used for:

- Moving through the menu
- Changing parameter value

OK button is used for:

- Selecting a menu from the display
- Confirming a chosen parameter value and leaving the current menu

NEXT button is used for going to the parameter of the current line.

PREVIOUS button is used for returning to the previous parameter of the current line.

MIXER ON/OFF button is used for starting/stopping the mixer.

TEMP ON/OFF button is used for starting/stopping the temperature regulator.

START/STOP button is used for starting/stopping the timer.

Display Contents

Power on the ENVIRO-GENIE (front rocker switch below the control panel). The display will briefly show the software revision.

The LCD display shows parameters in 4 lines.

Line 1:

- **MIXER MODE**
The ENVIRO-GENIE has 4 mixer modes.
Stirrer
Stirrer/Rev
Rotate/Rock
Off (Environmental Chamber - no mixing)
- **Rev. Time (sec)**
Displayed only when the mixer mode is **Stirrer/Rev** and programs the reverse stirring time in seconds.
- **Speed**
Stirrer or **Stirrer/Rev** mode: display shows XXX revolutions per minute.
Rotate/Rock mode: display shows XX/XX – the first set of values are for Rotate and the second set for Rock.
Off mode: display is blank.

Line 2 & 3:

- **TEMP**
MEAS: display shows XX.X °C – the current temperature in the environmental chamber.
SET: display shows XX.X °C – the desired temperature in the environmental chamber
- **ALARM:** display shows ±1°C thru ±9°C or NO. Alarm will sound if the chamber temperature deviates from your set temperature by the tolerance selected or will not sound if the alarm is disabled (NO).

Line 4:

- **TIMER:** display shows XX.XX.XX – the desired time in hours, minutes, seconds – counts down when cycle is started pressing ON/OFF button.
The Enviro-Genie has 4 **TIMER** modes:


OFF – expiration of time does not stop the mixer function or temperature regulation.

MIXER ONLY - expiration of time stops the mixer function.

TEMP ONLY – expiration of time stops the temperature regulation.

MIXER/TEMP – expiration of time stops both the mixer function and temperature regulation.

INCUBATOR GENIE

Button Function

UP and **DOWN** buttons in the **TEMP/°C** area are used for changing the temperature value.

UP and **DOWN** buttons in the **SPEED** area are used for changing the temperature value.

START/STOP button is used for starting/stopping rocking/rotating.

Display Contents

Power on the INCUBATOR GENIE (front rocker switch below the control panel). The display will briefly show the software revision.

- **MIXER MODE**

The INCUBATOR GENIE has 2 mixer modes.

Rotate/Rock

Off (Environmental Chamber - no mixing)

To activate **OFF** mixing mode, press the **SPEED DOWN** button until the display shows (--).

Mixing Modes Configuration

ENVIRO-GENIE & INCUBATOR GENIE

Rotate / Rock

The unit comes equipped with a magnetic platform that is used for rotating and rocking. Installing/removing the magnetic platform is easy but may take a few times to get used to. Try the following technique:

Install the magnetic platform by first unscrewing the lock shaft until the end of the lock shaft is flush with the right-hand sidewall of the inside of the chamber.

Grab the magnetic platform with the label side up, so that you can read the label. Hold the platform in the chamber and insert the hex shaft on the left side of the chamber into the hex hole in the magnetic platform. You may need to rotate the platform slightly to align the shaft with the hole.


Now hold up the right side of the platform, while depressing the blue lock release. Screw the lock shaft into the hole in the magnetic platform making sure the blue lock release is depressed. When the lock shaft is screwed all the way in, let go of the lock release.

To remove the platform, reverse the above procedure.


To select rotating mode:

1. Press the lock release.
2. Shift the magnetic platform to the right.
3. While shifting the platform right, release the lock release and line up the hex shaft to engage the lock shaft properly.
4. The platform will lock in place when it reaches the correct position.

**To select rocking mode:**

1. Press the lock release.
2. Shift the magnetic platform to the left.
3. While shifting the platform left, release the lock release and rotate the platform to line up the pin and the pin slot at the end of the platform. Shift the platform onto the pin.
4. The platform will lock in place when it reaches the correct position.

The versatile magnetic platform can be loaded in a variety of ways using both sides of the platform:
To rock open trays or racks, center on the metal rocker tray.

NOTE: The metal rocker tray must be used in rocking mode only. It will not clear the housing if rotated.

To rotate or rock bags, select mode of operation (Rotate or Rock); locate the bags in the desired position using 2 or 4 (depending on application) magnetic stainless steel mounting strips to hold them in place. Four mounting strips are supplied – more can be ordered as accessories. The platform capacity for bags can also be expanded.

Various size tubes can be placed on the magnetic platform in different orientations to produce different motions within the tubes. First select the size clip plate required. The large clips hold 28-30mm dia. (50ml) tubes, the medium clips hold 15-17mm dia. (15ml) tubes, and the small clips hold 10-13mm dia. (1.5-2.0ml) tubes. One clip plate of each size is supplied. The platform can hold up to eight plates. Push the tubes into the clips. Positioning the plates as shown on the right will produce a rotating or rocking motion in the tubes.

Positioning the plates as shown on the right will produce a rolling or rocking motion in the tubes.

Positioning the plates as shown on the right will produce a combination rotating/rolling motion ("Figure 8" motion) in the tubes. Adjust the angle to achieve the desired motion.

NOTE: When rotating tubes, rotate the magnetic platform by hand before starting, to check the clearance of the tubes to the chamber walls.

Environmental Chamber Only

To use the Environmental Chamber only, within the **MIXER MODE**, choose the setting **OFF**. All mixing will become disabled, and the unit becomes a temperature chamber only. A wire rack is included to increase the capacity of the chamber. An additional rack (sold separately) may be stacked on the first rack to further increase capacity.

Magnetic Stirring (ENVIRO-GENIE only)

ENVIRO-GENIE comes equipped with two magnetic stirring platforms, located at the bottom of the chamber. The stirrers are centered on the raised portions of the chamber bottom. The stirrers both always operate at the same settings. In general, the highest speeds can only be achieved with smaller stir bars with non-viscous liquids.

NOTE: If Rotating/Rocking Magnetic Platform is installed remove the platform prior to using magnetic stirrers.

Stir Reverse (ENVIRO-GENIE only)

Enabling “Stir Reverse” causes the stir bar to reverse direction at a rate set by the user. The stir bar is gently ramped down in speed to a stop, then ramped up to speed in the opposite direction. Time between reversing is adjustable from 1 to 60 seconds.

NOTE: If a very short stir reversing time is used in conjunction with a very fast stir speed, the commanded speed may not be attained. The system will not have time to ramp up to speed.

Operation Instructions

ENVIRO-GENIE

Mixer Operation

The mixer can be turned on by pressing the **MIXER ON/OFF** button or the **START/STOP** button (if parameters are set).

The mixer can be stopped by pressing the **MIXER ON/OFF** button or, if the time parameter expires, the unit will turn off automatically. While the mixer is not working, the **MIXER ON/OFF** button LED will not be lighted.

Temperature Regulator Operation

Temperature regulation can be turned ON/OFF by pressing the **TEMP ON/OFF** button or the **START/STOP** button. While the temperature regulation is working, the **TEMP ON/OFF** LED will be GREEN.

After the desired temperature has been reached, should the temperature deviate from the range set in the **ALARM** parameter, the **TEMP ON/OFF** LED will turn RED and the alarm will sound every 60 seconds. Once the desired temperature is achieved again, the LED will turn GREEN, and the alarm will stop.

Timer Operation

Timer can be switched on only by pressing the **START/STOP** button. While the timer is working, the **START/STOP** LED is GREEN. After the set time has expired, the **START/STOP** LED will blink, and an alarm will sound every 15 seconds. The time parameter will start to accumulate and represent the amount of time that has expired since the set time elapsed.

Pressing the **START/STOP** button will stop the time from accumulating, the **START/STOP** LED will turn off, and the time display will show the last value set in the time parameter.

Calibration Mode

This section applies **ONLY** to units with firmware version 1.20 and under. For units with version 1.21 or higher, see the “Calibration and PC Connection Manual” for calibration instructions.

The ENVIRO-GENIE is electronically factory calibrated and should not need calibration, nor is it recommended. However, should the need arise, the unit may be calibrated, as follows:

In order to calibrate, you will need to make an independent measurement of the chamber temperature. Make sure the thermometer is accurate! Generally speaking, glass thermometers found around the lab will be less accurate than your ENVIRO-GENIE. Good, certified glass or digital thermometers are available from laboratory supply dealers. A thermocouple wire may be passed through the door seal and will not affect performance.

Measure the temperature in the chamber via the lab thermometer, after chamber temperature has been allowed to settle for twenty minutes. Note the difference in temperature between ENVIRO-GENIE displayed temperature and your lab thermometer temperature (ex. ENVIRO-GENIE display reads 37.5°C, lab thermometer reads 37.0°C; difference is -0.5°C). This difference is the calibration value to be entered.

To access the calibration menu, turn the unit OFF.

Press and hold the **PREVIOUS** button while simultaneously turning the unit ON.

The parameters on the display show:

MEAS: XX.X – current temperature in the chamber read by regulatory temperature sensor.

SET: XX.X – calibration temperature setting in °C.

OFFSET: X.X – offset from the temperature setting.

The cursor on the display will blink at the value of the **SET** calibration point. Use the **UP** and **DOWN** buttons to adjust your desired **SET** point. After pressing the **OK** button, the temperature regulator will turn on and the cursor on the display will blink at the value of the **OFFSET** point. Wait until the **MEAS** and **SET** values are equal (approximately 20 minutes). Compare the temperature read by the unit temperature sensor, **MEAS**, and the temperature read by the lab thermometer. If the two values are different, adjust the **OFFSET**, using the **UP** and **DOWN** buttons, to make them equal. Wait 20 minutes. Repeat this procedure until the temperature sensor reading equals the lab thermometer reading.

Press the **OK** button. The display will move to the next menu.

Use the **UP** and **DOWN** buttons to choose:

SAVE&EXIT, **EXIT** or **RESTORE FACTORY**.

SAVE&EXIT saves the current value of the set temperature offset.

EXIT saves the old offset.

RESTORE FACTORY resets the set offset.

After pressing the OK button, the device restarts and goes to normal working mode.

If the **OFFSET** is 0°C, the display (top right) shows **FACTORY**. If the **OFFSET** is a value other than 0°C, the display shows **USER**.

RANGE OF OFFSET is -2°C to +2°C in 0.1°C increments.

RS-232 Interface

ENVIRO-GENIE comes equipped with a serial interface for connection to a PC or other serial device. This interface can be used to record machine information (including temperature) for situations that require verification. It can also be used to send commands remotely to the ENVIRO-GENIE, either manually or automatically from a program in order to create profiles.

The serial port is on the side of the unit and is a standard subminiature 9-pin male connector. Connection to a PC is via a straight through 9-pin female-female serial cable (sold separately).

- Baud Rate: 19,200
- Data Bits: 8
- Stop Bits: 1
- Parity: None

All communication is accomplished with ASCII Uppercase characters with a carriage return as the terminating character.

Control Commands

Name	Description
TIME	Mixing and temperature control time. Enter commanded time as follows: To enter 7 Hrs. 23 Min. type "Time 723" followed by the enter key.
TON	Turn on data output. Data is output once per second. To turn on data output, type the following: "TON" followed by the enter key. The format of the data out is as follows: Commanded Temperature, Time and Actual Temperature.
TOFF	Turn off data output. To turn off data output, type the following: "TOFF" followed by the enter key.
TEMP	Temperature command in degrees Celsius. Will accommodate tenths of degrees. To enter commanded temperature, type the following: "TEMP 250" followed by the enter key. This will return to 25.0 degrees Celsius on the front panel. "TEMP 251" followed by the enter key will return 25.1 degrees Celsius. The numerical value is the temperature value in degrees Celsius in the range from 4°C – 75°C.
SPEED	Speed command in RPM. To enter the commanded speed, type the following: "SPEED 500" followed by the enter key. The same ranges for Rock/Rotate (0-35 RPM) and Stirring (0-1000 RPM) entered on the keypad hold for RS-232 speed command.
MMODE	Mixing mode. To enter mixing mode, type the following: "MMODE 1-3" followed by the enter key. (1: No Mix mode. 2: Rock/Rotate mode. 3: Stir mode.)
TMODE	Timer mode. To enter timer mode, type the following: "TMODE 1-4" followed by the enter key. (1: Nothing off. 2: Mixing off. 3: Temperature off. 4: Mixing and temperature off.)

Name	Description
ALARM	Alarm mode. To enter alarm mode, type the following: "ALARM 1-3" followed by the enter key. (1: No temperature alarm. 2: $\pm 3^{\circ}\text{C}$ alarm. 3: $\pm 1^{\circ}\text{C}$ alarm.)
START	Start mixing/stirring. To start mixing or stirring, type the following: "START" followed by the enter key.
STOP	Stop mixing/stirring. To stop mixing or stirring, type the following: "STOP" followed by the enter key.
PELTON	To turn on heating/cooling control, type the following: "PELTON" followed by the enter key.
PELTOFF	To turn off heating/cooling control, type the following: "PELTOFF" followed by the enter key.
RESTOREE	Restore all factory-preset values.

WINDOWS Quick Start

1. Connect cable to ENVIRO-GENIE and PC.
2. Start Windows HyperTerminal (Start, Programs, Accessories, Communication, HyperTerminal).
3. Double Click on HyperTerminal.
4. Give the connection a name (i.e. "ENVIRO-GENIE").
5. Select the serial port you are connected to (usually COM1 or COM2).
6. Select 19,200, 8, None, 1.
7. You are now connected.

In order to have Windows scroll down after each line, check the box "append line feeds to incoming line ends."

Type commands to the unit in capital letters followed by the enter key.

In order to capture data in Windows HyperTerminal select "Transfer", "Capture Text" (Don't forget to issue the command "TON" to start data flowing from **ENVIRO-GENIE**).

INCUBATOR GENIE

Mixer Operation

The mixer can be turned on/off by pressing the **START/STOP** button.

Temperature Regulator Operation

When the unit is first turned **ON** (rocker switch below control panel), the heaters and fan will turn on and begin regulating toward the set temperature. During this time, the temperature display will blink green. The display will become solid green after achieving the desired temperature. Pressing the Up or Down Arrow buttons in the temperature control area allows you to change the desired temperature. While adjusting temperature, the display will be amber in color. Once the temperature has reached the set value, a change of $\pm 2^{\circ}\text{C}$ in the measured temperature will change the displayed value to solid red. At $\pm 3^{\circ}\text{C}$ deviation, the display will blink red. At $\pm 4^{\circ}\text{C}$ the display will blink red and the alarm will sound.

NOTE: This device is not capable of achieving a chamber temperature that is less than the ambient room temperature.

Calibration Mode

The INCUBATOR GENIE is electronically factory calibrated and should not need calibration, nor is it recommended. However, should the need arise, the unit may be calibrated, as follows:

In order to calibrate, you will need to make an independent measurement of the chamber temperature. Make sure the thermometer is accurate! Generally speaking, glass thermometers found around the lab will be less accurate than your INCUBATOR GENIE. Good, certified glass or digital thermometers are available from laboratory supply dealers. A thermocouple wire may be passed through the door seal and will not affect performance.

Measure the temperature in the chamber via the lab thermometer, after chamber temperature has been allowed to settle for twenty minutes. Note the difference in temperature between Incubator-Genie displayed temperature and your lab thermometer temperature (ex. Incubator-Genie display reads 37°C , lab thermometer reads 38°C ; difference is -1°C). This difference is the calibration value to be entered.

To access the calibration menu, simultaneously press and hold the **UP ARROW** button in the **TEMP/ $^{\circ}\text{C}$** area and the **START/STOP** button while turning the unit ON.

The parameters on the display show:

TEMPERATURE: The temperature blinks amber and begins adjusting to reach the last temperature setting. Once there, it changes to solid amber.

SPEED: 0 – The speed readout displays “0” in an amber color and can be adjusted ± 9 for each $^{\circ}\text{C}$ desired.

Press the **START/STOP** button to exit calibration mode.

Maintenance

The OHAUS MULTI-FUNCTION INCUBATOR GENIE should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. DO NOT use a cleaning agent or solvent on the front panel which is abrasive or harmful to plastics, nor one which is flammable. The unit and platform cannot and should not be autoclaved. The clip plates and mounting strips are made of magnetic stainless steel and can be cleaned with most cleaning agents or autoclaved. The metal rocker tray cannot and should not be autoclaved. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your Ohaus representative.

Technical Support Information

For service assistance or technical support, please visit our website at www.ohaus.com to locate and contact the Ohaus office nearest you.

Serial Number: _____

Date of Purchase: _____

Supplier: _____

Technical Data

Equipment Rating

- Indoor use only
- Altitude: 0 to 2000 m
- Environmental: 0°C to 38°C (32°F to 100°F), max 95% relative humidity, non-condensing
- Overvoltage category (Installation category): II

Compliance

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), (EU)2015/863, 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at www.ohaus.com/ce .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.ohaus.com/weee .
	EN 61326-1

Global Notice

Warning: This is a Class A product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

ISED Canada Compliance Statement:

CAN ICES-003(A) / NMB-003(A) FCC Notice

FCC Compliance Statement:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Enviro-Genie, e-MI3C... ESD and EFT C grade performance statement:

This product complies with IEC 61000-4-2:2008, Performance Criterion C. During severe Electrostatic Discharges (ESD) the equipment may experience a temporary loss of functionality or performance. If the device stops working or behaves abnormally, it will require user intervention, such as a restart, to restore normal operation.

ISO 9001 Registration

The management system governing the production of this product is ISO 9001 certified.

FCC Supplier Declaration of Conformity

Unintentional Radiator per 47CFR Part B

Trade Name: OHAUS CORPORATION

Model: ENVIRO-GENIE e-MI3C..., INCUBATOR GENIE e-MI3...

Party issuing Supplier's Declaration of Conformity

Troemner LLC

201 Wolf Drive

Thorofare, NJ, 08086

USA

Phone: +1 856 686 1600

Responsible Party – U.S. Contact Information

Ohaus Corporation

8 Campus Drive, Ste. 105

Parsippany, NJ 07054

United States

Phone: +1 973 377 9000

Web: www.ohaus.com

LIMITED WARRANTY

OHAUS products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period OHAUS will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to OHAUS.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than OHAUS. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by OHAUS Corporation. OHAUS Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact OHAUS or your local OHAUS dealer for further details.

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P/N: 31018510 A © 2026 Ohaus Corporation, all rights reserved