**Frontier™ 2000 Series
Centrifuges
Instruction Manual**

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LCD Display

The following picture shows the individual elements of the LCD-display.

Figure.3

Messages/logos of the display fields

“Prog”	“RPM/RCF”	“Time”	“Acc/Dec”
8 (Program)	8888 (RPM)	88:88 (Min:Sec)	88 (Accel Decel)
	88-88 (Rotor Number)	E 1 (Error Code)	
	8888. (RCF)		

Rotor Information Table

Rotor No. display	Order No.	Capacity	Compatible	
			FC2706	FC2516
03	30694903	Rotor, Angle, 12x1.5/2.0ml, ID, Sealable		●
04	30694904	Rotor, Angle, 10x5ml, ID, Sealable		●
08	30694908	Rotor, Angle, 24x1.5/2.0ml, ID, Sealable		●
09	30694909	Rotor, Angle, 6x10ml, ID, Sealable		●
44	30694844	Rotor, Angle, 48x1.5/2.0ml, ID, Sealable		●
38	30692838	Rotor, Angle, 6x15ml, ID	●	
07	30694907	Rotor, Angle, 12x15ml, ID	●	
11	30694811	Rotor, Angle, 6x50ml, ID	●	
14	30694814	Rotor, Angle, 6x100ml, ID	●	
65	30692865	Rotor, Swing out, 6x10ml, ID	●	

Table of Contents

1. SAFETY INFORMATION.....	1
1.1 Description and Intended Purpose.....	1
1.2 Brief description.....	1
1.3 Definition of Signal Warnings and Symbols	1
Safety Precautions	1
2.INSTALLATION.....	2
2.1 Unpacking.....	2
2.1.1 Delivery package.....	3
2.2 Selecting the Location	3
2.3 Installation	3
2.4 Safety precautions during operation and warranty	3
3.OPERATION	4
3.1 Mounting and loading rotor.....	4
3.1.1 Installation of rotors	4
3.1.2 Loading angle rotors.....	5
3.1.3 Loading swing out rotors.....	6
3.1.4 Loading and overloading of rotors	6
3.2 Lid control.....	6
3.2.1 Lid open	6
3.2.2 Lid lock	6
3.3 Parameter Settings.....	7
3.3.1 Rotor number setting	7
3.3.2 Rotation speed setting.....	7
3.3.3 Relative centrifugal force setting (RCF)	7
3.3.4 Centrifugation time setting	7
3.3.5 Acceleration rate setting	8
3.3.6 Deceleration rate setting.....	8
3.4 Programs.....	8
3.4.1 Saving programs	8
3.5 Starting and stopping the centrifuge	9
3.5.1 Starting the centrifuge	9
3.5.2 Stopping the centrifuge.....	9
3.6 Imbalance detection	9
4.MAINTENANCE	9
4.1 Maintenance and cleaning.....	9
4.1.1 General Care:.....	9
4.1.2 Cleaning and disinfection of the unit.....	10
4.1.3 Cleaning and disinfection of the rotor.....	10
4.1.4 Disinfection of aluminum rotors.....	10
4.1.5 Disinfection of PP-rotors	11
4.1.6 Glass breakage	11
4.1.7 Life time of rotors, buckets, accessories	11
5.TROUBLESHOOTING	11
5.1 Error message: Cause / Solution.....	11
5.2 Survey of possible error messages and their solutions	12
5.2.1 Lid release during power failure (Emergency Lid Release).....	12
5.2.2 Description of the error message system.....	12
6.RECEIPT OF CENTRIFUGES TO REPAIR.....	12
7.TRANSPORT AND STORAGE.....	12
7.1 Transport.....	12
7.2 Storage.....	12
7.3 Passing on the Device	13

8. TECHNICAL DATA.....	13
8.1 Specification.....	13
8.2 Drawings and dimensions.....	13
9. COMPLIANCE.....	14
10. APPENDIX.....	15
10.1 Table 1: Permissible net weight.....	16
10.2 Table 2: Max. speed and RCF-values for permissible rotors.....	16
10.3 Table 3: Acceleration and deceleration times	16

1. SAFETY INFORMATION

1.1 Description and Intended Purpose

Thank you for choosing this OHAUS product.

All symbols indicate safety instructions and points to potential dangerous situations. Please read the manual completely before using the Frontier™ 2000 Centrifuges to avoid incorrect operation.

Frontier™ 2000 Centrifuges were designed for the separation of materials or mixtures with different densities. OHAUS centrifuges are intended exclusively for indoor use and for use by qualified personnel.

1.2 Brief description

The model FC2516 is non-refrigerated micro high centrifuge.

The model FC2706 is non-refrigerated universal centrifuge.

All models are offered in two voltage variations, 230V and 120V.

The centrifuges can be used with swing-out and angle rotors.

The centrifuge is powered by a maintenance-free induction motor. Detailed technical data are in the "Technical data" section.

1.3 Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

Signal Words

WARNING	For a hazardous situation with medium risk, possibly resulting in injuries or death if not avoided.
CAUTION	For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or injuries if not avoided.
ATTENTION	For important information about the product. May lead to equipment damage if not avoided
NOTE	For useful information about the product

Warning Symbols

General Hazard

Explosion hazard

Electrical Shock Hazard

Biohazard

Safety Precautions

CAUTION: Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain all instructions for future reference.

- Before connecting power, verify that the AC adapter's input voltage range and plug type are compatible with the local AC mains power supply.
- Do not position the equipment such that it is difficult to reach the power connection.
- Only connect the power cord to a compatible grounded electrical outlet.
- Only use a power cord with a rating that exceeds the specifications on the equipment label.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- Operate the equipment only under ambient conditions specified in these instructions.
- This equipment is for indoor use only.
- Do not operate the equipment in wet, hazardous or unstable environments.

- Do not lean on the centrifuge.
- Do not stay within the clearance envelope longer than necessary for operational reasons.
- Do not allow liquids to enter the equipment.
- Do not place the equipment upside down on the platform.
- Disconnect the equipment from the power supply when cleaning.
- Service should only be performed by authorized personnel.
- Do not deposit any potentially hazardous materials within the clearance envelope.
- Do not unscrew the two halves of the housing.
- Only listed OHAUS original rotors and accessories for this instrument shall be used. Any other use or intended use is considered improper. OHAUS is not liable for damage resulting from improper use.
- During centrifugation, the presences of people are prohibited within 30 cm around the centrifuge.
- Under no circumstances open the lid of the centrifuge while the rotor is still running or rotating with a speed of $> 2\text{m/s}$.

WARNING: Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight. (explosion hazard due to spark formation, corrosion caused by the ingress of gases).

WARNING: Electrical shock hazards exist within the housing. The housing should only be opened by authorized and qualified personnel. Remove all power connections to the unit before opening.

WARNING: The centrifuge is not sealed. Use suitable protection measures when using the centrifuge for infectious and pathogenic samples. Follow appropriate safety precautions when handling these samples.

Danger and precautions

To protect people and environment the following precautions should be observed:

- During centrifugation, the presences of people are prohibited within 30 cm around the centrifuge according to the regulations of EN 61010-2-020.
- The equipment is not explosion-proof and must therefore not be operated in explosion-endangered areas or locations. Centrifugation of flammable, explosive, radioactive, or such substances, which chemically react with high energy, is strictly prohibited. The final decision on the risks associated with the use of such substances is the responsibility of the user of the centrifuge.
- Never spin toxic or pathogenic material without adequate safety precautions, i.e. centrifugation of buckets / tubes with missing or defective hermetic sealing is strictly prohibited. The user is obliged to perform appropriate disinfection procedures in case dangerous substances have contaminated the centrifuge and or its accessories. When centrifuging infectious substances, always pay attention to the general laboratory precautions. If necessary, contact your safety officer!
- It is prohibited to run the centrifuge with rotors other than listed for this unit.
- Under no circumstances open the lid of the centrifuge while the rotor is still running or rotating with a speed of $> 2\text{m/s}$.

2.INSTALLATION

2.1 Unpacking

Carefully remove your centrifuge and each of its components from the package. The included components vary depending on the centrifuge model (see table below). Save the packaging to ensure

safe storage and transport. The instruction manual must always be kept with the centrifuge!

Rotor(s) / Accessories will be packed separately.

WARNING: Lifting Hazard. Single person lift could cause injury. Use a mechanical lifting device or team lifting procedures when lifting or moving the equipment.

Please refer to section 8.3 for details about lifting the equipment out of the packaging.

2.1.1 Delivery package

Quantity	Description
1	Centrifuge FC2706/FC2516
1	Power Cable
1	Warranty Card
1	Instruction Manual / User Guide
1	Rotor Key

2.2 Selecting the Location

NOTE!

Avoid excessive vibrations, heat sources, air current, or rapid temperature changes.

- The centrifuge should be installed on an even, solid and level surface, if possible on a laboratory cabinet / table or some other solid vibration free surface.
- During centrifugation, the centrifuge must be placed in a way, that there is a minimum space of 30 cm on each side of the unit according to the standards EN 61010-2-020.
- Do not place the centrifuge next to a window or a heater, where it could be exposed to excessive heat, as the performance of the unit is based on an ambient temperature of 23°C.

2.3 Installation

Follow these steps:

- Check whether the power supply corresponds with the one specified on the manufacturer's rating label, which is located on the rear panel.
- For FC2706/FC2516, the power line should be protected by a 10 A rating circuit breaker (type K).
- In case of emergency, there must be an emergency switch off installed outside the room in order to disconnect the power supply from the unit.
- Connect the centrifuge to a grounded power receptacle.
- Connect the centrifuge with the mains. (The socket for the power cord must be easy to reach for disconnection)
- Turn the instrument on using the mains power switch.
- Open the lid by using the Door Open button.
- Remove the transport-securing device of the motor.

2.4 Safety precautions during operation and warranty

- Do not operate the centrifuge in case it's installed incorrectly.
- Do not lean on the centrifuge during operation.
- Do not stay within the 30 cm clearance envelope longer than necessary for operational reasons.

- Do not place any potentially hazardous materials within the 30 cm clearance envelope.
- Do not operate the centrifuge when disassembled (e.g. without housing).
- Do not run the centrifuge when mechanical or electrical components have been tampered with.
- Do not use accessories such as rotors and buckets, which are not exclusively approved by OHAUS Corporation, except commercially available centrifuge tubes made of glass or plastic.
- Do not spin extremely corrosive substances, as they may damage or weaken the materials.
- Do not operate the centrifuge with rotors or buckets, which show any signs of corrosion or mechanical damage.
- The manufacturer is responsible for safety and reliability of the centrifuge, only if:
 - 1) The unit is operated in accordance with this instruction manual.
 - 2) Modifications, repairs or other adjustments are performed by OHAUS authorized personnel and the electrical installation complies with the relevant electrical code.

NOTE!
Warranty

The centrifuge has been subjected to thorough testing and quality controls. In the unlikely case of any manufacturing faults occurring, the centrifuge and rotors are covered by warranty. This warranty becomes invalid in case of mishandling, damage and negligence and further in case of usage of inappropriate spare parts and / or accessories or unauthorized modification of the unit. Technical modification rights are reserved by the manufacturer in respect to technical improvement!

3. OPERATION

3.1 Mounting and loading rotor

3.1.1 Installation of rotors

Clean the drive shaft as well as the collet with a clean, grease-free piece of cloth. Place the rotor onto the drive shaft. (See figure below). Take care that the rotor is fully installed onto the motor shaft.

Motor shaft and chamber

Figure.4

Nut for Rotor

Tool for rotor
with nut

Nut for Rotor

Tool for rotor
with nut

Screw-on lid

Figure.5

Figure.6

Hold the rotor with one hand and secure the rotor to the shaft by turning the fixing nut clockwise. Tighten the fixing nut with enclosed rotor key (See figures 4-5)

We will provide a tool for nut rotor with centrifuge.

**ATTENTION!**

Check that the fixing screw is properly installed before each run. (See figure 5-6)

Do not operate the centrifuge with rotors or buckets which show any signs of corrosion or mechanical damage.

Do not operate with extremely corrosive substances, which could damage the rotor, buckets and materials.

In case of any questions, please contact the manufacturer!

3.1.2 Loading angle rotors

Rotors must be loaded symmetrically and with equal weight (See figure below). The adapter may only be loaded with the appropriate vessels. The weight differences between the filled vessels should be kept as low as possible. Therefore we recommend weighing them with a balance. This reduces the wear of the drive and the acoustic operating noise.

The maximum load per hole is stated on each rotor.

Figure.7 WRONG

Figure.8 CORRECT (6 tubes)

3.1.3 Loading swing out rotors

Loading of the buckets / vessels must be made in accordance to the figure below.

It is allowed to operate e.g. a 4-place-rotor with 2 loaded buckets only. But the loaded buckets must be opposite to each other. Make sure that the unloaded buckets also be put inside the rotor (see below).

In principle swing out rotors may not be taken into operation until all buckets or racks are put into the rotor. The bolts at the rotor must be greased with the "High TEF oil". The sample tubes have to be filled evenly by eye and put into the drillings or tube racks. The weight difference of the loaded buckets should not exceed approx. 1.0 g.

ATTENTION!

Swing out rotors may be taken in operation only if all locations are filled in with either four buckets or four carriers – do not mix buckets and carriers together!!

ATTENTION!

Do not operate the centrifuge with rotors or buckets which show any signs of corrosion or mechanical damage.

Do not operate with extremely corrosive substances, which could damage the rotor and buckets. In case of any questions, please contact the manufacturer!

Figure.9 WRONG

Figure.10 CORRECT

3.1.4 Loading and overloading of rotors

All approved rotors are listed with their maximum speed and maximum filling weight in "**table 2 permissible net weight**" (See APPENDIX).

The maximum load permitted for a rotor, which is determined by the manufacturer, as well as the maximum speed allowed for this rotor (See label on rotor), must not be exceeded. The liquids the rotors are loaded with should have a maximum homogeneous density of 1.2 g/ml or less when the rotor is running at maximum speed.

In order to spin liquids with a higher density, the speed has to be reduced according to the following formula:

$$\text{Reduced speed } n_{\text{red}} = \sqrt{\frac{1.2}{\text{higher density}}} \times \text{max. speed } (n_{\text{max}}) \text{ of the rotor}$$

Example:

$$n_{\text{red}} = \sqrt{\frac{1.2}{1.7}} \times 4.000 = 3.360 \text{ rpm}$$

If in case of any questions, please contact the manufacturer!

To determine the relative centrifugal force(RCF/G-force) for a specific adapter, you can calculate per DIN 58 970 using the attached formula:

$$\text{RCF} = 1.117862 \cdot 10^{-5} \cdot n^2 \cdot r_{\text{max}}$$

n: revolutions per minute (RPM)

r_{max}: max centrifuging radius in cm by using the bottom of tubes

3.2 Lid control

3.2.1 Lid open

Press the  button. After the mechanical lock is unlocked, the lid can be opened. Note: The top cover will automatically pop open with a sound when unlocked.

3.2.2 Lid lock

Press down slightly harder on the lid. The lid will close after the mechanical lock connects to the lock buckle.

3.3 Parameter Settings

3.3.1 Rotor number setting

Press the "SET" button until the RPM/RCF window displays a flashing "NO-38". Use the "▲" or "▼" buttons to set the correct rotor number. The rotor number is the last two digits of the rotor material serial number. E.g., the serial number of rotor 6×15ml is "30692838". So, the rotor number is "38". The setting is shown in the figure below:

Figure.11

3.3.2 Rotation speed setting

Press the "SET" button until the RPM/RCF window displays a flashing "8888". Use the "▲" or "▼" buttons to adjust the values. Pressing the "▲" or "▼" buttons once will adjust the rate by 5 rpm, while pressing down on either button for approximately 3 seconds will adjust the rate by 100 rpm. The figure below shows the setting at 4000 rpm.

Figure.12

3.3.3 Relative centrifugal force setting (RCF)

Press the "SET" button until the RPM/RCF window displays a flashing "8888.". Use the "▲" or "▼" buttons to adjust the values. Pressing the "▲" or "▼" buttons once will adjust the rate by 5×g, while pressing down on either button for approximately 3 seconds will adjust the rate by 100×g. The figure below shows the setting at 3500×g.

Figure.13

3.3.4 Centrifugation time setting

Minutes setting: Press the "SET" button until a flashing "88:" is displayed in the Time window. Use the "▲" or "▼" buttons to adjust the values by 1 minute each time you press. The range of values is 0 min-99 min.

Seconds setting: Press the "SET" button until a flashing ":88" is displayed in the Time window. Use the "▲" or "▼" buttons to adjust the values by 1 second each time you press. The range of values is 1s-59s.

The centrifuge timer will begin counting down after the rotation speed reaches the configured speed.

The figure below shows the configuration at 10 min 20s.


Figure.14

3.3.5 Acceleration rate setting

Press the "SET" button until a flashing "8" is displayed in the Acc window. Use the "▲" or "▼" buttons to adjust the value by 1 gear each time you press. The range of acceleration rate is 0- 9 gears. The acceleration rate of gear 0 is the slowest while the rate of gear 9 is the fastest. The higher the value, the faster the acceleration rate.

3.3.6 Deceleration rate setting

Press the "SET" button until a flashing "8" is displayed in the Dec window. Use the "▲" or "▼" buttons to adjust the value by 1 gear each time you press. The range of deceleration rate is 0 - 9 gears. The deceleration rate of gear 0 is the slowest while the rate of gear 9 is the fastest. The higher the value, the faster the deceleration rate. The figure below shows the acceleration rate at the gear 3 and the deceleration rate at gear 6.


Figure.15

3.4 Programs

3.4.1 Saving programs

You can save at most 9 sets of A-L3 programs, including operating parameters and the rotors used. Of these, A/B/C/D/E/F/G are single-step operating programs, H1-H2 are dual-step programs, and L1-L2-L3 are triple-step programs.

After installing the rotor, press the "Prog" button to select the program number that needs to be set. After setting the rotor number and operating parameters, wait for approximately 3 seconds and the set parameters will be automatically saved to the current program number.

3.4.2 Calling programs

Press the "Prog" button to select a pre-set program number (as shown at M1 in the Figure). The saved operating parameters will be displayed automatically.


Figure.16

3.5 Starting and stopping the centrifuge

3.5.1 Starting the centrifuge

Press the "Start" button and the instrument begins to operate. The green indicator light on the button will illuminate, the rotating speed will gradually rise to the set value and the countdown timer will begin to work. If it is necessary to adjust the rotation speed during the centrifuging process, press the "Set" button until the rotation speed shown in the RPM/RCF window is flashing. Then use the "▲" or "▼" buttons to adjust the value. After that, the instrument will operate according to the newly set rotation speed.

Warning: While the machine is in operation, the operator shall not lean on the instrument. Non-staff personnel shall not linger in the safe space.

3.5.2 Stopping the centrifuge

(1) As the instrument keeps running, the countdown timer gradually turns to zero and the red indicator light on the stop button will illuminate. The instrument will begin decelerating from the operational rotation until it fully stops. Upon that, a buzzer will sound, indicating that the instrument has stopped operating and the lid can be opened.

(2) If the instrument has to be stopped during operation, press the  button. A red signal light on the button will illuminate. The instrument will fully stop as stated above.

Note: It is prohibited to manually unlock the instrument if the rotor speed has not decelerated to a safe range. Manually forcing the rotor to fully stop is strictly prohibited.

3.6 Imbalance detection

In case of the rotor not being equally loaded, the drive will turn off during acceleration. The rotor decelerates to stand still.

When in the display "Time" the word "E9" (M9), the weight difference of the samples is too large. Distribute the weight evenly.


Figure.17

4. MAINTENANCE

4.1 Maintenance and cleaning

4.1.1 General Care:

Maintenance of the centrifuge is confined to keeping the rotor, the rotor chamber and the rotor accessories clean as well as to regularly lubricating the rotor insert bolts of a swing out rotor (if available).

The most suitable lubricant is the High TEF oil.

Lubricants containing molycote and graphite are not allowed.

Please pay special attention to anodized aluminum parts. Breakage of rotors can be caused even by slight damage.

In case of rotors, buckets or tube racks getting in touch with corrosive substances the concerned spots have to be cleaned carefully.

Corrosive substances are for instance: alkalis, alkaline soap solutions, alkaline amines, concentrated acids, solutions containing heavy metals, water-free chlorinated solvents and saline solutions. e.g. salt water, phenol, halogenated hydrocarbons.

Cleaning – units, rotors, accessories:

- Turn the device off and disconnect it from the power supply before you begin any cleaning or disinfecting. Do not pour liquids into the housing interior.
- Do not spray disinfectant on the device.
- Thorough cleaning not only has its purpose in hygiene but also in avoiding corrosion due to pollution.
- In order to avoid damaging anodized parts such as rotors, reduction plates etc., only pH-neutral Detergents with a pH-value of 6-8 may be used for cleaning. Alkaline cleaning agents (pH-value > 8) must not be used.
- After cleaning, please ensure all parts are dried thoroughly, either by hand or in a hot-air cabinet (max. Temperature + 50°C).
- It is necessary to coat anodized aluminum parts with anti-corrosion oil regularly in order to increase their life- spans and reduce corrosion predisposition.
- Due to humidity or not hermetically sealed samples, condensate may be formed. The condensate has to be removed from the rotor chamber with a soft cloth regularly.

The maintenance procedure has to be repeated every 10 to 15 runs, or at least once a week.

- Connect the unit to the power supply, after the equipment is completely dry.
- Do not carry out disinfection with UV-, beta- and gamma-rays or other high energy radiation.
- Metal rotors can be autoclaved.
- Rotor lid and adapters can also be autoclaved (max. 121°C, 20 min).
- The tube racks are made of PP and **cannot** be autoclaved at 134°C.

4.1.2 Cleaning and disinfection of the unit

1. Open the lid before you turn off the unit. Disconnect it from the power supply.
2. Open the rotor nut by turning the rotor key counter clockwise.
3. Remove the rotor.
4. For cleaning and disinfection of the unit and the rotor chamber use the above mentioned cleaner.
5. Clean all accessible areas of the device and its accessories, including the power cord with a damp cloth.
6. Wash the rubber seals and rotor chamber thoroughly with water.
7. Rub the dry rubber seals with glycerol or talc to prevent these to becoming brittle. Other components of the unit, e.g. the lid lock, motor shaft and rotor must not be greased.
8. Dry the motor shaft with a soft, dry and lint-free cloth.
9. Control the unit and accessories for damage.

Make sure that the centrifuge is turned off the unit and disconnect the unit from the power supply. Then remove adherent dust from the ventilation slots in the centrifuge by using a soft brush. Do this at least every six months.

WARNING:

Electric Shock Hazard. Disconnect the equipment from the power supply before cleaning. Make sure that no liquid enters the interior of the equipment.

4.1.3 Cleaning and disinfection of the rotor

1. Clean and disinfect the rotors, rotor lids and adapters with the above mentioned cleaner.
2. Use a bottle brush to clean and disinfect the rotor bores.
3. Rinse the rotors, rotor lid and adapter with clear water. Particularly the drillings of angle rotors.
4. For drying of the rotors and accessories set them on a towel. Place the angle rotors with bores down.
5. Dry the rotor cone with a soft, dry and lint-free cloth and look for damage. Do not grease the rotor cone.

4.1.4 Disinfection of aluminum rotors

In case of infectious material spilling into the centrifuge, the rotor and rotor chamber have to be disinfected

directly after the run. Rotors may be autoclaved at a maximum temperature of 121°C.

4.1.5 Disinfection of PP-rotors

Autoclaving

The recommended time for autoclaving: 15 – 20 min at 121°C (1 bar)

ATTENTION:

The sterilization time of 20 min. must not be exceeded. Repeated sterilization will cause reduction of the mechanical resistance of the plastic material

Before autoclaving the PP-rotor and adapter must be thoroughly cleaned to avoid the burning in of dirty residues. You can disregard the consequences of some chemical residues to plastic materials at ambient temperatures. But at the high temperatures during autoclaving those residues may corrode and destroy the plastic. The objects must be thoroughly rinsed with distilled water after the cleaning but before the autoclaving. Residues of any cleaning liquids may cause fissures, whitening and stains.

Gas sterilization

Adapters, bottles and rotors may be gas sterilized with Ethylenoxyd. Make sure to air out the items after the sterilization and before using them again.

ATTENTION:

Because the temperature may rise during the sterilization, rotors, adapters and bottles must not be closed and must be totally unscrewed.

Chemical sterilization

Bottles, adapters and rotors may be treated with the usual liquid disinfectants.

ATTENTION:

Before applying any other cleaning or decontamination method than recommended by the manufacturer, contact the manufacturer to ensure that it will not damage the unit or the rotor.

4.1.6 Glass breakage

With high g-values, the rate of glass tube breakage increases. Glass splinters have to be removed immediately from rotor, buckets, adapters and the rotor chamber itself. Fine glass splinters will scratch and therefore damage the protective surface coating of a rotor. If glass splinters remain in the rotor chamber, fine metal dust will build up due to air circulation. This very fine, black metal dust will significantly pollute the rotor chamber, the rotor, the buckets and the samples.

If necessary, replace the adapters, tubes and accessories to avoid further damage. Check the rotor bores regularly for residues and damage.

ATTENTION:

Please check the relevant specifications of the tubes centrifuges with the manufacturer.

4.1.7 Life time of rotors, buckets, accessories

Rotors and rotor lid made of aluminum or stainless steel, have an operating time of max. 7 years from first use. Transparent rotor lids and caps made of PC or PP as well as rotors, tube racks and adapters of PP have a maximum operating time up to 3 years from first use.

Condition for the operating time: Proper use damage-free condition, recommended care.

5.TROUBLESHOOTING

5.1 Error message: Cause / Solution

The error messages are listed to help localize possible errors faster.

The diagnosing referred to in this chapter may not always be the case, as they are only theoretically occurring errors and solutions.

Please keep us informed about any kind of error occurring, which is not listed in this chapter. Only through your information are we able to improve this operation manual.

Many thanks in advance for your support.

5.2 Survey of possible error messages and their solutions

5.2.1 Lid release during power failure (Emergency Lid Release)

The OPEN button suddenly fails due to shutdown as a result of power failure or due to shutdown as a result of the burnout of fuse of the instrument. In order to protect the centrifuge sample, a pull ring is set at the bottom of the instrument. Please pull the pull ring downward, the door lock will be opened, and the door cover can be pushed open. However, the following rule-breaking operations are strictly prohibited.

(1) In the case the revolutions per minute of the instrument is not reduced to 100r/min during speed reduction and shutdown, it is strictly prohibited to pull the pull ring at the bottom of the instrument to forcibly open the door lock and manually force the rotor to stop operation, as which will cause personal injury.

(2) It is strictly prohibited to manually turn off the power switch and then turn on the power to press the OPEN button to open the lid and force the rotor to stop by hand at the point when the program of the instrument considers that the revolutions per minute is 0.

5.2.2 Description of the error message system

When operates the instrument improperly or a failure occurs, an error code will appear in the "Time" window. See the following table for information about error codes.

Error code	Description
E1	The door is not closed properly or a failure occurs to the door lock
E2	Control system failure
E3	Overvoltage fault
E4	Overcurrent fault
E7	Speed sensor fault
E8	Overspeed fault
E9	Unbalance fault

6.RECEIPT OF CENTRIFUGES TO REPAIR

Health risk from contaminated equipment, rotors and accessories.

In case of returning the centrifuge for repairing to the manufacturer, please notice the following:

The centrifuge must be decontaminated and cleaned before the shipment for the protection of persons, environment and material.

Decontamination certificate at goods return delivery (See APPENDIX). We reserve the right to not accept contaminated centrifuges.

Further on all costs occurred for the cleaning and disinfection of the units will go to the debit of the customer's account.

7.TRANSPORT AND STORAGE

7.1 Transport

Before transporting, take out the rotor.

Only transport the unit in the original packaging.

Install the transport protection material to secure the motor shaft, when transporting over longer distances.

	Air temperature	rel. humidity	Air pressure
General transportation	-25 to 60 °C	10 to 75 %	30 to 106 kPa

7.2 Storage

During storage of the centrifuge the following environmental conditions must be observed:

	Air temperature	rel. humidity	Air pressure
in transport packaging	-25 to 55 °C	10 to 75 %	70 to 106 kPa

7.3 Passing on the Device

When passing the equipment on to third parties, please make sure to also include this instruction manual.

8. TECHNICAL DATA

8.1 Specification

Equipment Ratings:

Indoor use only

Altitude: 2000m

Operating temperature: 15 to 35°C

Humidity: Maximum relative humidity 80% for temperatures up to 31 °C decreasing linearly to 50% relative humidity at 35°C.

Electrical supply: 230V ~, 4A, 50/60Hz.


Voltage fluctuations: Mains supply voltage fluctuations up to $\pm 10\%$ of the nominal voltage.


Overvoltage category (Installation category): II

Pollution degree: 2

Model	FC2706	FC2516
Speed Range	300 rpm - 6000 rpm; 5 rpm/set	300 rpm -16500 rpm; 5 rpm/set
Maximum RCF	3660 x g; 5 x g/set	19080 x g; 5 x g/set
Maximum Capacity(Rotor)	6x100 ml	48x1.5/2 ml, 6x10 ml
Temperature range(N/A)	Air cool	
Running Time	1 sec to 99 min 59 sec	
Noise level (depending on the rotor)	$\leq 63 \pm 2$ dB(A)	
Allowable density at maximum speed	1.2 g/ml	
Mains power connection AC	230 V ~ 50/60 Hz	
Voltage fluctuation	$\pm 10\%$	
Current consumption	3 A	4 A
Power consumption	300 W	400 W
Net Weight (without rotor)	21 KG 46.3 lb	15.5 KG 34.2 lb
Shipping Dimensions (W x D x H)	53 x 40 x 35.5 cm	47 x 36 x 35.5 cm
Shipping Weight (without rotor)	24 KG	18 KG

8.2 Drawings and dimensions

Model	W (mm / in.)	L (mm / in.)	H (mm / in.)
FC2706	320 / 12.6	440 / 17.3	250 / 9.8
FC2516	280 / 11.0	380 / 15.0	245 / 9.7

9.COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

Mark	Standard
	This product complies with the applicable harmonized standards of EU Directives 2011/65/EU (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at www.ohaus.com/ce .
	This product complies with the applicable statutory standards of the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012, UK Electromagnetic Compatibility Regulations 2016 and Electrical Equipment (Safety) Regulations 2016. The UK Declaration of Conformity is available online at www.ohaus.com/uk-declarations .
	This product complies with the EU Directive 2012/19/EU (WEEE). Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. For disposal instructions in Europe, refer to www.ohaus.com/weee .
	EN 61326-1

ISO 9001 Registration
The management system governing the production of this product is ISO 9001 certified.

10.APPENDIX

TABLE 1: PERMISSIBLE NET WEIGHT

TABLE 2: MAX. SPEED AND RCF-VALUES FOR PERMISSIBLE ROTORS

TABLE 3: ACCELERATION AND DECELERATION TIMES

10.1 Table 1: Permissible net weight

Rotor number	Rotor Description	Permissible net weight	Max Speed (×1000 rpm)	
			FC2706	FC2516
30692838	Rotor Angle 6×15ml	6 x 22.5g	6.0	
30694907	Rotor Angle 12×15ml	12 x 22.5g	5.0	
30694811	Rotor Angle 6×50ml	6 x 75g	5.0	
30694814	Rotor Angle 6×100ml	6 x 150g	5.0	
30692865	Rotor Swing Out 6×10ml	6 x 15g	4.0	
30694903	Rotor Angle 12x1.5/2ml	12 x 3g		16.5
30694904	Rotor Angle 10x5ml	10 x 7.5g		16.0
30694908	Rotor Angle 24x1.5/2.0ml	24 x 3g		14.0
30694909	Rotor Angle 6×10ml	6 x 15g		13.0
30694844	Rotor Angle 48x1.5/2ml	48 x 3g		12.0

10.2 Table 2: Max. speed and RCF-values for permissible rotors

Rotor number	Rotor Description	Permissible net weight	Model	Max Speed	
				RPM	RCF
30692838	Rotor Angle 6×15ml	6 x 22.5g	FC2706	6000	3660×g
30694907	Rotor Angle 12×15ml	12 x 22.5g	FC2706	5000	3080×g
30694811	Rotor Angle 6×50ml	6 x 75g	FC2706	5000	2850×g
30694814	Rotor Angle 6×100ml	6 x 150g	FC2706	5000	3130×g
30692865	Rotor Swing Out 6×10ml	6 x 15g	FC2706	4000	2200×g
30694903	Rotor Angle 12x1.5/2ml	12 x 3g	FC2516	16500	19080×g
30694904	Rotor Angle 10x5ml	10 x 7.5g	FC2516	16000	17882×g
30694908	Rotor Angle 24x1.5/2.0ml	24 x 3g	FC2516	14000	17956×g
30694909	Rotor Angle 6×10ml	6 x 15g	FC2516	13000	14196×g
30694844	Rotor Angle 48x1.5/2ml	48 x 3g	FC2516	12000	14079×g

10.3 Table 3: Acceleration and deceleration times

Rotor number	Rotor Description	Permissible net weight	Model	Acceleration* Time in sec		Deceleration* Time in sec	
				level 0	level 9	level 0	level 9
30692838	Rotor Angle 6×15ml	6 x 22.5g	FC2706	68	24	39	14
30694907	Rotor Angle 12×15ml	12 x 22.5g	FC2706	60	22	84	15
30694811	Rotor Angle 6×50ml	6 x 75g	FC2706	59	22	180	15
30694814	Rotor Angle 6×100ml	6 x 150g	FC2706	59	20	180	15
30692865	Rotor Swing Out 6×10ml	6 x 15g	FC2706	48	21	38	13
30694903	Rotor Angle 12x1.5/2ml	12 x 3g	FC2516	110	45	159	23
30694904	Rotor Angle 10x5ml	10 x 7.5g	FC2516	109	44	184	23
30694908	Rotor Angle 24x1.5/2.0ml	24 x 3g	FC2516	96	40	280	23
30694909	Rotor Angle 6×10ml	6 x 15g	FC2516	88	37	355	25
30694844	Rotor Angle 48x1.5/2ml	48 x 3g	FC2516	82	38	492	32

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P/N 30602522B © 2022 Ohaus Corporation, all rights reserved.
Ohaus Corporation