**Frontier™ Centrifuge  
FC5706  
Instruction Manual**

**Centrifuga Frontier™  
FC5706  
Manual de instrucciones**

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Front and rear view of the centrifuge FC5706


Figure.1

|                     |                    |
|---------------------|--------------------|
| 1 Centrifuge Lid    | 2 Rotor Window     |
| 3 Display           | 4 Function Label   |
| 5 Main Power Switch | 6 Power Connection |
| 7 Emergency Release |                    |

Function Label

Function Label for FC5706


Figure.2

|                                                        |                                          |
|--------------------------------------------------------|------------------------------------------|
| 1 LCD Display                                          | 2 Start centrifugation                   |
| 3 Stop centrifugation / setup                          | 4 RPM/RCF model and select               |
| 5 Acceleration/Deceleration intensity model and select | 6 Time setup model                       |
| 7 Release lid                                          | 8 Adjusting knob/Dial: Change the number |

## LCD Display

The following picture shows the individual elements of the LCD-display.


Figure.5

### Display fields:

- A-1 Display field – "RPM/RCF"
- A-2 Display field – "Acc/Dec" "Service"
- A-4 Display field – "Temp"

### Messages/logos of the display fields

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## TABLE OF CONTENTS

|                                                                      |    |
|----------------------------------------------------------------------|----|
| 1. INTRODUCTION .....                                                | 1  |
| 1.1 Description .....                                                | 1  |
| 1.2 Features .....                                                   | 1  |
| 1.3 Definition of Signal Warnings and Symbols .....                  | 1  |
| 1.4 Safety Precautions .....                                         | 2  |
| 1.4.1 User .....                                                     | 2  |
| 1.4.2 Rotor and accessories .....                                    | 2  |
| 1.4.3 Measures for your protection .....                             | 2  |
| 1.4.4 Exclude the following environmental influences .....           | 2  |
| 1.4.5 Measures for operational safety .....                          | 3  |
| 1.4.6 Danger and precautions .....                                   | 3  |
| 1.4.7 Abbreviations used in this manual .....                        | 3  |
| 2. INSTALLATION .....                                                | 4  |
| 2.1 Unpacking .....                                                  | 4  |
| 2.1.1 Delivery package for Frontier™ FC5706 .....                    | 4  |
| 2.2 Selecting the Location .....                                     | 4  |
| 2.3 Installation .....                                               | 4  |
| 2.4 Safety precautions during operation .....                        | 4  |
| 3. OPERATION .....                                                   | 5  |
| 3.1 Mounting and loading rotor .....                                 | 5  |
| 3.1.1 Installation of rotors .....                                   | 5  |
| 3.1.2 Loading angle rotors .....                                     | 5  |
| 3.1.3 Loading swing out rotors .....                                 | 5  |
| 3.1.4 Loading and overloading of rotors .....                        | 6  |
| 3.1.5 Removing the rotor .....                                       | 6  |
| 3.2 Power switch .....                                               | 6  |
| 3.3 Lid control .....                                                | 7  |
| 3.3.1 Lid open .....                                                 | 7  |
| 3.3.2 Lid lock .....                                                 | 7  |
| 3.4 Preselection .....                                               | 8  |
| 3.4.1 Preselection of speed / .....                                  | 8  |
| 3.4.2 Preselection of running time .....                             | 8  |
| 3.4.3 Preselection of brake intensity and acceleration .....         | 9  |
| 3.5 Starting and stopping the centrifuge .....                       | 10 |
| 3.5.1 Starting the centrifuge .....                                  | 10 |
| 3.5.2 The "STOP" key .....                                           | 10 |
| 3.6 Imbalance detection .....                                        | 11 |
| 4. SETTING .....                                                     | 11 |
| 4.1 Change the type of rotor .....                                   | 11 |
| 4.2 Access to mode "Operating Data" .....                            | 12 |
| 4.3 Call up operating data .....                                     | 12 |
| 5. MAINTENANCE .....                                                 | 13 |
| 5.1 Maintenance and cleaning .....                                   | 13 |
| 5.1.1 General Care. ....                                             | 13 |
| 5.1.2 Cleaning and disinfection of the unit .....                    | 14 |
| 5.1.3 Cleaning and disinfection of the rotor .....                   | 14 |
| 5.1.4 Disinfection of aluminum rotors .....                          | 14 |
| 5.1.5 Disinfection of PP-rotors .....                                | 14 |
| 5.1.6 Glass breakage .....                                           | 14 |
| 5.2 Lifetime of rotors, buckets, accessories .....                   | 15 |
| 6. TROUBLESHOOTING .....                                             | 15 |
| 6.1 Error message: Cause / Solution .....                            | 15 |
| 6.2 Survey of possible error messages and their solutions .....      | 15 |
| 6.2.1 Lid release during power failure (Emergency Lid Release) ..... | 15 |
| 6.2.2 Description of the error message system .....                  | 15 |
| 7. RECEIPT OF CENTRIFUGES TO REPAIR .....                            | 16 |
| 8. TRANSPORT, STORAGE AND DISPOSAL .....                             | 16 |
| 8.1 Transport .....                                                  | 16 |
| 8.2 Storage .....                                                    | 16 |
| 9. TECHNICAL DATA .....                                              | 17 |
| 9.1 Specifications .....                                             | 17 |

|      |                                                                 |    |
|------|-----------------------------------------------------------------|----|
| 9.2  | Drawings and dimensions.....                                    | 18 |
| 10.  | ORDER INFORMATION.....                                          | 18 |
| 10.1 | Rotor .....                                                     | 18 |
| 11.  | COMPLIANCE.....                                                 | 19 |
| 12.  | APPENDIX .....                                                  | 20 |
| 12.1 | Table 1: Permissible net weight.....                            | 20 |
| 12.2 | Table 2: Max. speed and RCF-values for permissible rotors ..... | 20 |
| 12.3 | Table 3: Acceleration and deceleration times.....               | 20 |
| 12.4 | Table 4: Error messages .....                                   | 21 |
| 12.5 | Table 5: Redemption form / Decontamination certificate.....     | 22 |

## 1. INTRODUCTION

### 1.1 Description

Thank you for choosing this OHAUS product.

All symbols indicate safety instructions and points to potential dangerous situations. Please read the manual completely before using the Frontier™ FC5706 to avoid incorrect operation.

Frontier™ 5706 centrifuge was designed for the separation of materials or mixtures with different density.

### 1.2 Features

The Frontier™ FC5706 centrifuge offers many practical features such as:

- Distinct control panel
  - Simple one-handed operation
  - Set and view speed in both rpm and g-force
- Control all settings with ergonomic adjustment knob
  - Accommodate all tube sizes
  - 1.5ml, 2.0ml, 5ml, 7ml, 15ml, 16ml, 30ml, 50ml
- No tool needed to install or uninstall rotors

### 1.3 Definition of Signal Warnings and Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results.

The degree of danger is a part of a safety note and distinguishes the possible results of non-observance from each other.

#### Signal Words

|                  |                                                                                                                                            |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>DANGER</b>    | Will lead to severe injuries or death if not avoided.                                                                                      |
| <b>WARNING</b>   | For a hazardous situation with medium risk, possibly resulting in injuries or death if not avoided.                                        |
| <b>CAUTION</b>   | For a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or injuries if not avoided. |
| <b>ATTENTION</b> | For important information about the product. May lead to equipment damage if not avoided                                                   |
| <b>NOTE</b>      | For useful information about the product                                                                                                   |

#### Warning Symbols

General Hazard

Electrical Shock Hazard

Alternating Current

Biohazard

Explosion

Crushing

## Warning and information signs on the surface of centrifuge

### Warning

Four carrier must be used at all times on four place swing out rotors or damage will occur to the centrifuge. Such damage will not be covered under the product warranty.

Four carrier must be used at all times on four place swing out rotors or damage will occur to the centrifuge. Such damage will not be covered under the product warranty.

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**Attention! Check the fastening of the rotor nut before each run.**

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TAKE OFF MAINS PLUG before opening  
the housing or the emergency release!

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**Take off mains plug before opening the housing or the emergency release.**

## 1.4 Safety Precautions

### 1.4.1 User

OHAUS centrifuges are intended exclusively for indoor use and for use by qualified personnel. This device may only be operated by trained specialist staff. They must have carefully read the operating manual and be familiar with the function of the device.

### 1.4.2 Rotor and accessories

Only OHAUS original rotors and accessories shall be used. Any other use or intended use is considered improper. OHAUS is not liable for damage resulting from improper use.

#### CAUTION:

Read all safety warnings before installing, making connections, or servicing this equipment. Failure to comply with these warnings could result in personal injury and/or property damage. Retain instructions for future reference.

### 1.4.3 Measures for your protection

**WARNING:** Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight. (Explosion hazard due to spark formation, corrosion caused by the ingress of gases)

**WARNING:** When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules.

**WARNING:** The centrifuge is not sealed. Use suitable protection measures when using the centrifuge for infectious and pathogenic samples. Follow appropriate safety precautions when handling these samples.

### 1.4.4 Exclude the following environmental influences

- Powerful vibrations
- Direct sunlight
- Atmospheric humidity greater than 80%
- Corrosive gases present
- Temperatures below 2 °C and above 35 °C
- Powerful electric or magnetic fields:

#### WARNING:

Electrical shock hazards exist within the housing. The housing should only be opened by authorized and qualified personnel.

Remove all power connections to the unit before opening.

**1.4.5 Measures for operational safety**

- Do not unscrew the two halves of the housing.
- Dry off any liquid spills immediately! The instrument is not watertight.
- Verify that the equipment's input voltage range and plug type are compatible with the local power supply.
- Only connect the power cord to a properly grounded power receptacle.
- Only use a power cord with a rating that exceeds the specifications on the equipment label.
- Do not position the equipment such that it is difficult to disconnect the power cord from the power receptacle.
- Make sure that the power cord does not pose a potential obstacle or tripping hazard.
- The equipment is for indoor use only. Use the equipment only in dry locations.
- Use only approved accessories.
- Operate the equipment only under ambient conditions specified in these instructions.
- Disconnect the equipment from the power supply when cleaning.
- Do not operate the equipment in hazardous or unstable environments.
- Service should only be performed by authorized personnel.

**1.4.6 Danger and precautions**

To protect people and environment the following precautions should be observed:

- During centrifugation, the presences of people are prohibited within 30 cm around the centrifuge according to the regulations of EN 61010-2-020.
- FC5706 is not explosion-proof and must therefore not be operated in explosion-endangered areas or locations. Centrifugation of flammable, explosive, radioactive, or such substances, which chemically react with high energy, is strictly prohibited. The final decision on the risks associated with the use of such substances is the responsibility of the user of the centrifuge.
- Never spin toxic or pathogenic material without adequate safety precautions, i.e. centrifugation of buckets / tubes with missing or defective hermetic sealing is strictly prohibited. The user is obliged to perform appropriate disinfection procedures in case dangerous substances have contaminated the centrifuge and or its accessories. When centrifuging infectious substances, always pay attention to the general laboratory precautions. If necessary, contact your safety officer!
- It is prohibited to run the centrifuge with rotors other than listed for this unit.
- Under no circumstances open the lid of the centrifuge while the rotor is still running or rotating with a speed of > 2m/s.

**1.4.7 Abbreviations used in this manual**

| Symbol/Abbreviations | Unit                       | Description                   |
|----------------------|----------------------------|-------------------------------|
| RPM                  | $[\text{min}^{-1}]$<br>rpm | revolutions per minute        |
| RCF                  | [x g]                      | relative centrifugal force    |
| PCR                  |                            | PCR Polymerase chain reaction |
| PP                   | -                          | Polypropylene                 |
| PC                   | -                          | Polycarbonate                 |
| accel                | -                          | acceleration                  |
| decel                | -                          | deceleration                  |

## 2. INSTALLATION

### 2.1 Unpacking

Carefully remove your centrifuge and each of its components from the package. The included components vary depending on the centrifuge model (see table below). Save the packaging to ensure safe storage and transport. The instruction manual must always be kept with the centrifuge!

**Rotor(s) / Accessories will be packed separate.**

#### 2.1.1 Delivery package for Frontier™ FC5706

| Quantity | Description        |
|----------|--------------------|
| 1        | Centrifuge FC5706  |
| 1        | Power Cable        |
| 1        | Warranty Card      |
| 1        | Instruction Manual |

### 2.2 Selecting the Location

#### Attention!

Avoid excessive vibrations, heat sources, air current, or rapid temperature changes.

- The centrifuge should be installed on an even, solid and level surface, if possible on a laboratory cabinet / table or some other solid vibration free surface.
- During centrifugation, the centrifuge must be placed in a way, that there is a minimum space of 30 cm on each side of the unit according to the standards EN 61010-2-020.
- Do not place the centrifuge next to a window or a heater, where it could be exposed to excessive heat, as the performance of the unit is based on an ambient temperature of 23°C.

### 2.3 Installation

Follow these steps:

- Check whether the power supply corresponds with the one specified on the manufacturer's rating label, which is located on the rear panel.
- The line voltage circuit breaker is max. 10 A (type K) slow release for commonly used instruments.
- In case of emergency, there must be an emergency switch off installed outside the room in order to disconnect the power supply from the unit.
- Connect the centrifuge to a grounded power receptacle.
- Turn the instrument on using the mains power switch.
- Open the lid by using the Door Open button.
- Remove the transport securing device of the motor.

### 2.4 Safety precautions during operation

- Do not operate the centrifuge in case it is not installed correctly.
- Do not lean on the centrifuge during operation.
- Do not stay within the 30 cm clearance envelope longer than necessary for operational reasons.
- Do not place any potentially hazardous materials within the 30 cm clearance envelope.
- Do not operate the centrifuge when disassembled (e.g. without housing).
- Do not run the centrifuge when mechanical or electrical components have been tampered with.
- Do not use accessories such as rotors and buckets, which are not exclusively approved by OHAUS Corporation, except commercially available centrifuge tubes made of glass or plastic.
- Do not spin extremely corrosive substances, as they may damage or weaken the materials.
- Do not operate the centrifuge with rotors or buckets, which show any signs of corrosion or mechanical damage.

**The manufacturer is responsible for the safety and reliability of the centrifuge, only if:**

- The unit is operated in accordance with this instruction manual.
- Modifications, repairs or other adjustments are performed by authorized personnel and the electrical installation complies with the relevant electrical code.

### 3. OPERATION

#### 3.1 Mounting and loading rotor

##### 3.1.1 Installation of rotors

Clean the drive shaft as well as the collet with a clean, grease-free piece of cloth. Place the rotor onto the drive shaft. (See figure below) Take care that the rotor is fully installed onto the motor shaft.

Figure.10


Figure 11

Figure 12

Hold the rotor with one hand and secure the rotor to the shaft by turning the fixing screw clockwise. (See figure 12)

##### ATTENTION:

Check that the fixing screw is properly installed before each run.

Do not operate the centrifuge with rotors or buckets which show any signs of corrosion or mechanical damage.

Do not operate with extremely corrosive substances, which could damage the rotor, buckets and materials.

In case of any questions, please contact the manufacturer!

##### 3.1.2 Loading angle rotors

Rotors must be loaded symmetrically and with equal weight (See figure below). The adapter may only be loaded with the appropriate vessels. The weight differences between the filled vessels should be kept as low as possible. Therefore we recommend weighing them with a balance. This reduces the wear of the drive and the acoustic operating noise.

On each rotor, the maximum load per hole is stated. (It is only allowed to operate e.g. a 12-place-rotor with 2 ,4 or 8 loaded tubes. But the loaded borings must be opposite each other).

Figure.13-1 WRONG

Figure.13-2 CORRECT (6 tubes)

##### 3.1.3 Loading swing out rotors


Loading of the buckets / vessels must be made in accordance with figure 16. It is only allowed to operate e.g. a 4-place-rotor with 2 loaded buckets, with the loaded buckets located opposite to each other and the unloaded buckets put inside the rotor (See figure below).

In principle, swing out rotors may not be taken into operation until all buckets or racks are put into the rotor.

The sample tubes have to be filled evenly by eye and put into the drillings or tube racks. The weight difference of the loaded buckets should not exceed 1 g.

**ATTENTION!**

Swing out rotors may be taken into operation only if all locations are filled in with either four buckets or four carriers – do not mix buckets and carriers up!

Figure. 16-1: **WRONG**Figure. 16-1: **CORRECT****ATTENTION!**

Do not operate the centrifuge with rotors or buckets which show any signs of corrosion or mechanical damage.

Do not operate with extremely corrosive substances, which could damage the rotor and buckets.

In case of any questions, please contact the manufacturer!

**3.1.4 Loading and overloading of rotors**

All approved rotors are listed with their maximum speed and maximum filling weight in **"table 2 permissible net weight"** (See APPENDIX).

The maximum load permitted for a rotor, which is determined by the manufacturer, as well as the maximum speed allowed for this rotor (See label on rotor), must not be exceeded. The liquids the rotors are loaded with should have a maximum homogeneous density of 1.2 g/ml or less when the rotor is running at maximum speed.

In order to spin liquids with a higher density, the speed has to be reduced according to the following formula:

$$\text{Reduced speed } n_{\text{red}} = \sqrt{\frac{1,2}{\text{high} \diamond \diamond}} \times \text{max. speed } (n_{\text{max}}) \text{ of the rotor}$$

Example:

$$\frac{\diamond \diamond \diamond \diamond \diamond \diamond \diamond \diamond}{y}$$

$$n_{\text{red}} = \sqrt{\frac{1,2}{1,7}} \times 4,000 = 3,360 \text{ rpm}$$

In case of any questions, please contact the manufacturer!

**3.1.5 Removing the rotor**

Untighten the rotor fixing nut completely (screw over the stiff point) and lift the rotor vertically out of the centrifuge. (See figure 12)

**3.2 Power switch**

The power switch is located on the bottom left side of the unit (see figure 17).

Figure. 17: Power switch

### 3.3 Lid control

#### 3.3.1 Lid open

After the run, when closing the lid of the centrifuge, in the display "**RPM | RCF**" (A-1) the word "**close**" (M1) appears (refer to figure 18 below).

If there is a rotor in the centrifuge, additional the word "**rotor**" (M3) appears, as well as the code number of the respective rotor which is in the centrifuge system "**77**" (M4). If there is no rotor in the centrifuge, the word "**rotor**" (M3) flashes and additional the word "**no**" (M4) appears. By pressing the key "**Door Open**" (7) you can release the lid of the centrifuge. As soon as the electromagnetic lid is completely released, the word "**open**" (M2) appears. Now you can open the lid of the centrifuge.


Figure. 18

During the run you can call up the rotor type at any time by pressing the key "**Door Open**" (7).

#### 3.3.2 Lid lock

The lid should only be put down slightly. An electromagnetic lid lock closes the lid, at the same time the word "**open**" (M2) disappears (refer to figure 18).

As a sign that the centrifuge is ready for starting, in the display "**RPM | RCF**" (A-1) the word "**close**" (M1) appears. Simultaneously the word "**rotor**" (M3) is displayed, as well as the code number of the rotor, which is in the centrifuge system, "**no 71**" (M4). With that, all rotor specific data, like e. g. max. speed, acceleration etc., are adopted.

#### ATTENTION:

Before closing the lid please check if the rotor is tighten, and all 6 buckets have put in the swing out rotor.

### 3.4 Preselection

#### 3.4.1 Preselection of speed / RCF-value

This pre-selection is activated through the key **"RPM | RCF"** (4) (refer to figure 19 below). By pressing the key once the word **"rpm"** (M5) flashes. By pressing the key twice the pre-selection of the centrifugal forces can be selected. Then the flashing word **"rcf"** (M6) appears. You can set the desired values with the adjusting knob (1). In the display (A-1) the regulated value is shown permanently, before, during and after the run.


Figure. 19

The speed is adjustable between 200 rpm and maximum revolution of the centrifuge resp. the maximum permissible revolution of the pre-selected rotor.

It is the same as the pre-selection of the RCF-value. The setting range is between the centrifugal force at 200 rpm and the maximum permissible centrifugal force of the rotor.

The maximum speed of the FC5706 is 6000 rpm resp. 4427 x g.

See **"Table 3: max. speed and RCF-values for permissible rotor"** (See APPENDIX). All important values are listed there.

#### ATTENTION:

Please also check the maximum permissible revolutions of your test tubes with the manufacturer.

#### 3.4.2 Preselection of running time

The running time can be pre-selected in three different ranges from 10 seconds up to 99 hours 59 minutes.

1. Range from 10 seconds up to 59 minutes 50 seconds in steps of 10 seconds
2. Range from 1 hour up to 99 hours 59 minutes in steps of 1 minute
3. The continuous run **"cont"**, can be interrupted by the key **"Stop"** (10).

The running time can be pre-selected with the lid open or closed.

To activate the setting of the running time press the key **"Time"** (6).

In the display **"Time"** (A-3) the indication **"m : s"** or **"h : m"** flashes, depending on the previous setting.

Set the desired value by using the adjusting knob (1). After exceeding 59 min 50 sec the indication changes automatically into **"h : m"**. After exceeding 99 hours 59 min the word **"cont"** appears in the display **"Time"** (A-3). That continuous run can only be interrupted by pressing the key **"Stop"** (10). The time starts counting down as soon as the set speed is reached.

The display always shows the remaining running time. (See figure 20)

All with number marked passages refer to figure 20


Figure. 20

### 3.4.3 Preselection of brake intensity and acceleration

This function is activated through the key **"Accel/Decel"** (5) (refer to figure 21).

By pressing the key once the word **"accel"** (M7) flashes in the display **"Acc/Dec"** (A-2). The desired acceleration can be pre-selected by the adjusting knob (1). The value 0 is equivalent to the lowest and the value 9 to the highest acceleration.

By pressing the key **"Accel/Decel"** (5) twice, in the display **"Acc/Dec"** (A-2), the word **"decel"** (M8) is shown. Now the desired brake intensity can be pre-selected by the adjusting knob (1). The value 9 is equivalent to the shortest and the value 0 to longest possible brake time.

See **"table 4: acceleration and deceleration times"** (APPENDIX). There the acceleration and deceleration times for the acceleration and deceleration stages 0 to 9 for permissible rotors are shown.


Figure. 21

## 3.5 Starting and stopping the centrifuge

### 3.5.1 Starting the centrifuge

To start the centrifuge, make sure that the lid is closed and then press the key **"Start"** (9), (refer to figure 22). With the key **"Start"** (9) you can start runs with manually pre-selected parameters. When the respective pre-selected running time has ended then the centrifuge will stop automatically or you can interrupt each run with the key **"Stop"** (10).


Figure. 22

### 3.5.2 The "STOP" key

The run can be interrupted at any time by pressing the **"Stop"** key (10) (See figure 23). After pressing the key the centrifuge decelerates with the respective pre-selected intensity until it reaches standstill.


Figure. 23

### 3.6 Imbalance detection

In case of the rotor not being equally loaded, the drive will turn off during acceleration. The rotor will decelerate to a standstill.

When in the display **"Time"** (A-3) the word **"error"** (M11) together with the number **"01"** appear, the weight difference of the samples is too large. Weigh the samples more exactly and start again.

Load the rotor as described in chapter 3.1.2 and 3.1.3.

When inside the display **"Time"** (A-3) the word **"error"** together with the number **"02"** (See figure 31) appear, it could be due to the imbalance switch being defective.


Figure. 24

## 4. SETTING

### 4.1 Change the type of rotor

Before the first operation and after each rotor change, you have to set the rotor type. You find each rotor type in the printed order number on the rotor.

Example:

Angle rotor order number: 30130877

Rotor Type on the display = 77

Turn on the centrifuge and open the lid. Now simultaneously press the keys **"Door Open"**(7) and **"Stop"**(10). In the display **"RPM | RCF"** the old rotor type no **"77"** appears. With the potentiometer you can now set the rotor type. To store the new setting please press the **"Start"** (9) key. Inside the display, **"Store"** appears as confirmation.

Now all important rotor parameters for the centrifuge are stored.

#### Attention!

The set rotor type must always be the same as the actual rotor type used; otherwise the equipment might be damaged.

The rotor type can be checked during the run by pressing the key **"Door Open"**(7).


Figure. 25

## 4.2 Access to mode "Operating Data"

When using the centrifuge, the following parameters can be set:

- Number of starts
- Operating hours of centrifuge
- Operating hours of motor
- Software-version
- Error list
- Function of the imbalance switch
- Operation of keyboard
- Display tests

While the centrifuge is turned off, press simultaneously the keys **"Time"**(6) and **"Door Open"** (7) and turn on the main switch of the centrifuge. Now release both keys and as a result a display test is executed for approx. 5 seconds. All indicators will appear at the same time (See figure 26).


Figure. 26

### ATTENTION:

After the settings have been stored, change back to normal mode again by switching off the centrifuge for a short while!

## 4.3 Call up operating data

### ATTENTION:

This should only be performed by advance user or service engineer.

In the mode **"Basic Adjustments"** you can call up the operating data of the centrifuge. Please proceed as described under point 4.1.2 to enter this program mode. Press the key **"Accel/Decel"** (5). In the display **"Acc/Dec"** (A-2) the word **"Service"** flashes.

With the adjusting knob (1) various information can be called up:

A= previous starts of the centrifuge

H= previous operating hours

S= software version

r= converter software

E= list of previous error messages

h= running time of the motor

The list of the last 99 error messages can accessed by pressing the key **"RPM | RCF"** (4). Scroll through the list by using the adjusting knob (1). The respective error codes appear in the display **"RPM | RCF"** (A-1). Please refer to **"Table 5: error messages"** (see APPENDIX).


Figure. 27

## 5. MAINTENANCE

### 5.1 Maintenance and cleaning

#### 5.1.1 General Care.

Maintenance of the centrifuge is confined to keeping the rotor, the rotor chamber and the rotor accessories clean as well as to regularly lubricating the rotor insert bolts of a swing out rotor (if available).

The most suitable lubricant is the offered OHAUS High TEF oil – Order no.: 30130896.

Lubricants containing molycot and graphite are not allowed.

Please pay special attention to anodized aluminum parts. Breakage of rotors can be caused even by slight damage.


In case of rotors, buckets or tube racks getting in touch with corrosive substances the concerned spots have to be cleaned carefully.

Corrosive substances are for instance: alkalis, alkaline soap solutions, alkaline amines, concentrated acids, solutions containing heavy metals, water-free chlorinated solvents, saline solutions, e.g. salt water, phenol, halogenated hydrocarbons.

#### Cleaning – units, rotors and accessories

- Turn the device off and disconnect it from the power supply before you begin any cleaning or disinfecting. Do not pour liquids into the housing interior.
- Do not spray disinfectant on the device.
- Thorough cleaning not only has its purpose in hygiene but also in avoiding corrosion due to pollution.
- In order to avoid damaging anodized parts such as rotors, reduction plates etc., only pH-neutral Detergents with a pH-value of 6-8 may be used for cleaning. Alkaline cleaning agents (pH-value > 8) must not be used.
- After cleaning, please ensure all parts are dried thoroughly, either by hand or in a hot-air cabinet (max. Temperature + 50°C).
- It is necessary to coat anodized aluminum parts with anti-corrosion oil regularly in order to increase their life-spans and reduce corrosion predisposition.
- Due to humidity or not hermetically sealed samples, condensate may be formed. The condensate has to be removed from the rotor chamber with a soft cloth regularly.

The maintenance procedure has to be repeated every 10 to 15 runs, or at least once a week.

- Connect the unit to the power supply, after the equipment is completely dry.
- Do not carry out disinfection with UV-, beta- and gamma-rays or other high energy radiation.
- Metal rotors can be autoclaved.
- Rotor lid and adapters can also be autoclaved (max. 121°C, 20 min).
- The tube racks are made of PP and cannot be autoclaved at 134°C.

### 5.1.2 Cleaning and disinfection of the unit

1. Open the lid before you turn off the unit. Disconnect it from the power supply.
2. Open the rotor nut by turning the rotor key counter clockwise.
3. Remove the rotor.
4. Cleaning and disinfect the unit and the rotor chamber using the above mentioned cleaner.
5. Clean all accessible areas of the device and its accessories, including the power cord with a damp cloth.
6. Wash the rubber seals and rotor chamber thoroughly with water.
7. Rub the dry rubber seals with glycerol or talc to prevent these from becoming brittle. Other components of the unit, e.g. the lid lock, motor shaft and rotor must not be greased.
8. Dry the motor shaft with a soft, dry and lint-free cloth.
9. Check the equipment and accessories for damage.

Make sure that the centrifuge is turned off the unit and disconnect the unit from the power supply. Then remove adherent dust from the ventilation slots in the centrifuge by using a soft brush. Do this at least every six months.

### 5.1.3 Cleaning and disinfection of the rotor

1. Clean and disinfect the rotors, rotor lids and adapters with the above mentioned cleaner.
2. Use a bottle brush to clean and disinfect the rotor bores.
3. Rinse the rotors, rotor lid and adapter with clear water. Particular the drillings of angle rotors.
4. For drying of the rotors and accessories set them on a towel. Place the angle rotors with bores down, to dry them to.
5. Dry the rotor cone with a soft, dry and lint-free cloth and look for damage. Do not grease the rotor cone.

### 5.1.4 Disinfection of aluminum rotors

In case of infectious material spilling into the centrifuge, the rotor and rotor chamber have to be disinfected right after the run. Rotors may be autoclaved at a maximum temperature of 121 °C.

### 5.1.5 Disinfection of PP-rotors

#### Autoclaving

The recommended time for autoclaving: 15 – 20 min at 121 °C (1 bar)

#### **ATTENTION:**

The sterilization time of 20 min. must not be exceeded. Repeated sterilization will cause reduction of the mechanical resistance of the plastic material  
Before autoclaving the PP-rotor and adapter must be thoroughly cleaned to avoid the burning in of dirty residues.

You can disregard the consequences of some chemical residues to plastic materials at ambient temperatures. But at the high temperatures during autoclaving those residues may corrode and destroy the plastic. The objects must be thoroughly rinsed with distilled water after the cleaning but before the autoclaving. Residues of any cleaning liquids may cause fissures, whitening and stains.

#### Gas sterilization

Adapters, bottles and rotors may be gas sterilized with Ethylenoxyd. Make sure to air out the items after the sterilization and before using them again.

#### **ATTENTION:**

Because the temperature may rise during the sterilization, rotors, adapters and bottles must not be closed and must be totally unscrewed.

#### Chemical sterilization

Bottles, adapters and rotors may be treated with the usual liquid disinfectants.

#### **ATTENTION:**

Before applying any other cleaning or decontamination method than recommended by the manufacturer, contact the manufacturer to ensure that it will not damage the unit or the rotor.

### 5.1.6 Glass breakage

With high g-values, the rate of glass tube breakage increases. Glass splinters have to be removed immediately from rotor, buckets, adapters and the rotor chamber itself. Fine glass splinters will scratch and therefore damage the protective surface coating of a rotor. If glass splinters remain in the rotor chamber, fine metal dust will build up due to air circulation. This very fine, black metal dust will seriously pollute the rotor chamber, the rotor, the buckets and the samples.

If necessary, replace the adapters, tubes and accessories to avoid further damages. Check the rotor bores regularly for residues and damage.

**ATTENTION:**

Please check the relevant specifications of the tubes with the manufacturer.

## 5.2 Lifetime of rotors, buckets, accessories

Rotors and rotor lid made of aluminum or stainless steel, have an operating time of max. 7 years from first use. Transparent rotor lids and caps made of PC or PP as well as rotors, tube racks and adapters of PP have a maximum operating time up to 3 years from first use.

Condition for the operating time: Proper use, damage-free condition and recommended care.

## 6. TROUBLESHOOTING

### 6.1 Error message: Cause / Solution

The error messages are listed to help localize possible errors faster.

The cause and solution referred to in this chapter may not always be the case, as they are only theoretically occurring errors and solutions.

Please keep us informed about any kind of error occurring, which is not listed in this chapter. Only through your information are we able to improve and this operation manual.

Many thanks in advance for your support.

### 6.2 Survey of possible error messages and their solutions

#### 6.2.1 Lid release during power failure (Emergency Lid Release)

In case of power failure or malfunction, the lid of the centrifuge can be opened manually in order to protect your samples.

Please proceed as follows:

1. Switch the centrifuge off and unplug the power cord, wait until the rotor stands still (this may take several minutes).
2. At the left side of the centrifuge housing there is a plastic stopper (see figure 28). Remove this stopper, fastened to it there is a string which is connected to the electronic lid lock.
3. If you pull the string slightly the lid will open.

**ATTENTION:**

- Don't put your hands in the rotor chamber as long as the rotor is still spinning!
- Push the plastic stopper back in the unit again, to continue working.

Figure. 28

#### 6.2.2 Description of the error message system

The error message **"error"** (M11) is shown in the **"Time"** (A-3) display (See figure 29). Detailed information about possible error messages are in **"table 5: error messages"** (See Appendix).


Figure. 29

## 7. RECEIPT OF CENTRIFUGES TO REPAIR

Health risk from contaminated equipment, rotors and accessories.

In case of returning the centrifuge for repairing to the manufacturer, please notice the following:

The centrifuge must be decontaminated and cleaned before the shipment for the protection of persons, environment and material.

A decontamination certificate needs to be filled in and enclosed with any returned goods. (See APPENDIX)

We reserve the right to not accept contaminated centrifuges.

Furthermore, all costs occurred for the cleaning and disinfection of the units will go to the debit of the customer's account.

## 8. TRANSPORT, STORAGE AND DISPOSAL

### 8.1 Transport

Before transport, take out the rotor.

Only transport the unit in the original packaging.

Use a transport aid for transporting over longer distances to fix the motor shaft.

|                        | Air temperature | rel. humidity | Air pressure  |
|------------------------|-----------------|---------------|---------------|
| General transportation | -25 to 60 °C    | 10 to 75 %    | 30 to 106 kPa |

### 8.2 Storage

During storage of the centrifuge the following environmental conditions must be observed:

|                        | Air temperature | rel. Humidity | Air pressure  |
|------------------------|-----------------|---------------|---------------|
| in transport packaging | -25 to 55 °C    | 10 to 75 %    | 70 to 106 kPa |


## 9. TECHNICAL DATA

### 9.1 Specifications

|                                                 |                                                                                                                  |                  |
|-------------------------------------------------|------------------------------------------------------------------------------------------------------------------|------------------|
| Model                                           | FC5706                                                                                                           |                  |
| Speed Range                                     | 200 rpm -6000 rpm;50 rpm/set                                                                                     |                  |
| Maximum RCF                                     | 4427 x g;10 x g/set                                                                                              |                  |
| Maximum Capacity(Rotor)                         | 6 x 50 ml                                                                                                        |                  |
| Temperature range                               | Air cool                                                                                                         |                  |
| Running Time                                    | 10 sec to 99 hr 99 min 59 sec or continuous                                                                      |                  |
| Noise level (depending on the rotor)            | ≤ 60 +2 dB(A)                                                                                                    |                  |
| Allowable density at maximum speed              | 1.2 g/ml                                                                                                         |                  |
| Allowable kinetic energy                        | 2427 Nm                                                                                                          |                  |
| Mains power connection AC                       | 230 V ~ 50/60 Hz                                                                                                 | 120 V ~ 50/60 Hz |
| Voltage fluctuation                             | ± 10 %                                                                                                           | ± 10 %           |
| Current consumption                             | 0.55 A                                                                                                           | 1.1 A            |
| Power consumption                               | 100 W                                                                                                            | 100 W            |
| Dimensions (W × D × H)                          | 291 x 392 x 266 mm<br>11.5 x 15.4 x 10.5 in                                                                      |                  |
| Net Weight (without rotor)                      | 10.5 kg<br>23 lb                                                                                                 |                  |
| Shipping Dimensions (W × D × H)                 | 400 x 545 x 410 mm<br>15.7 x 21.5 x 16.1 in                                                                      |                  |
| Shipping Weight (without rotor)                 | 15 kg<br>33 lb                                                                                                   |                  |
| Ambient conditions (EN/IEC 61010-1)             |                                                                                                                  |                  |
| Environment                                     | for indoor use only                                                                                              |                  |
| Altitude                                        | Use up to an altitude of 2000 m                                                                                  |                  |
| Ambient temperature                             | 2℃ up to 35 ℃                                                                                                    |                  |
| Max. relative humidity                          | Max. relative humidity 80 % for temperatures up to 31℃, decreasing linearly to 50 % relative humidity up to 35℃. |                  |
| Overvoltage category (IEC 60364-4-443)          | II                                                                                                               |                  |
| Degree of contamination                         | 2                                                                                                                |                  |
| Class of protection                             | I                                                                                                                |                  |
| Not suitable for use in hazardous environments. |                                                                                                                  |                  |
| EMC                                             | EN/IEC 61326-1 Class B emissions, Basic immunity<br>FCC Class B emissions                                        |                  |

## 9.2 Drawings and dimensions

Dimensions for FC5706

## 10. ORDER INFORMATION

### 10.1 Rotor

| Rotor/Adapter   | Description                         | Units    |
|-----------------|-------------------------------------|----------|
| <b>30130877</b> | <b>Rotor Angle 12x15ml RB/FA</b>    | <b>1</b> |
| 30211180        | Adapter 1x9-10ml D16mm Vac/Sar 2/pk | 6        |
| 30130889        | Adapter 1x7ml D13mm RB/Vac/Sar 2/pk | 6        |
| 30130890        | Adapter 1x5ml D13mm RB/Vac/Sar 2/pk | 6        |
| 30130886        | Adapter 1x1.5/2.0ml D11mm 6/pk      | 2        |
| <b>30130878</b> | <b>Rotor Angle 6x50ml RB/FA</b>     | <b>1</b> |
| 30130891        | Adapter 1x30ml D26mm 2/pk           | 3        |
| 30130892        | Adapter 1x16ml D18mm 2/pk           | 3        |
| 30130893        | Adapter 1x15ml D17mm 2/pk           | 3        |
| 30130894        | Adapter 1x15ml D17mm RB/FA 2/pk     | 3        |
| <b>30130880</b> | <b>Rotor Swing out 6x5ml RB</b>     | <b>1</b> |

Note: Package is dependent on country and might vary. Please check the item number with the local OHAUS office before you order it.

## 11. COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

| Marking                                                                           | Standard                                                                                                                                                                                                               |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | This product conforms to the EMC Directive 2014/30/EU and the Low Voltage Directive 2014/35/EU. The complete Declaration of Conformity is available online at <a href="http://www.ohaus.com/ce">www.ohaus.com/ce</a> . |

|                                                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <p><b>Disposal</b></p> <p>In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.</p> <p>Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.</p> <p>If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.</p> <p>Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.</p> <p>For disposal instructions in Europe, refer to <a href="http://www.ohaus.com/weee">www.ohaus.com/weee</a>.</p> <p>Thank you for your contribution to environmental protection.</p> |
|-----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### FCC Note

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## 12. APPENDIX

TABLE 1: PERMISSIBLE NET WEIGHT

TABLE 2: MAX. SPEED AND RCF-VALUES FOR PERMISSIBLE ROTORS

TABLE 3: ACCELERATION AND DECELERATION TIMES

TABLE 4: ERROR MESSAGES

TABLE 5: REDEMPTION FORM / DECONTAMINATION CERTIFICATE

### 12.1 Table 1: Permissible net weight

| Order No. | Capacity                  | Max. speed | Permissible net weight |
|-----------|---------------------------|------------|------------------------|
| 30130877  | Rotor Angle 12x15ml RB/FA | 6000 rpm   | 300 g                  |
| 30130878  | Rotor Angle 6x50ml RB/FA  | 6000 rpm   | 432 g                  |
| 30130880  | Rotor Swing out 6x5ml RB  | 4000 rpm   | 240 g                  |

### 12.2 Table 2: Max. speed and RCF-values for permissible rotors

| Order No. | Capacity                  | Max. speed | RCF      |
|-----------|---------------------------|------------|----------|
| 30130877  | Rotor Angle 12x15ml RB/FA | 6000 rpm   | 4427 x g |
| 30130878  | Rotor Angle 6x50ml RB/FA  | 6000 rpm   | 4427 x g |
| 30130880  | Rotor Swing out 6x5ml RB  | 4000 rpm   | 1878 x g |

### 12.3 Table 3: Acceleration and deceleration times

| Order No.  | Capacity                  | Acceleration values                                               |         | Deceleration values                                               |         |
|------------|---------------------------|-------------------------------------------------------------------|---------|-------------------------------------------------------------------|---------|
|            |                           | Stage 0                                                           | Stage 9 | Stage 0                                                           | Stage 9 |
| 30130877   | Rotor Angle 12x15ml RB/FA | 95                                                                | 33      | 257                                                               | 51      |
| 30130878   | Rotor Angle 6x50ml RB/FA  | 91                                                                | 41      | 274                                                               | 52      |
| 30130880   | Rotor Swing out 6x5ml RB  | 35                                                                | 8       | 25                                                                | 7       |
| in seconds |                           |                                                                   |         |                                                                   |         |
|            |                           | Acceleration time<br>from 0 min <sup>-1</sup> -> U <sub>max</sub> |         | Deceleration time<br>from U <sub>max</sub> -> 0 min <sup>-1</sup> |         |

**12.4 Table 4: Error messages**

| Error-no.: | Description                                       |
|------------|---------------------------------------------------|
| 1          | Imbalance                                         |
| 2          | Imbalance sensor is defective                     |
| 14         | Leap of speed is too big between two measurements |
| 30         | Motor is blocked or defective                     |
| 33         | Open lid during the rotor is running              |
| 55         | Overspeed                                         |
| 60         | Undervoltage in the intermediate circuit          |
| 70         | sticking relay                                    |

## 12.5 Table 6: Redemption form / Decontamination certificate

Enclose this form with all returns of equipment and assemblies!

The complete declaration about the decontamination is a prerequisite for further processing of the return. If no corresponding explanation is enclosed, we carry out decontamination at the expense of the customer.

Surname; last name: \_\_\_\_\_

Organization / company: \_\_\_\_\_

Street: \_\_\_\_\_

ZIP CODE: \_\_\_\_\_ City: \_\_\_\_\_

Telephone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-Mail: \_\_\_\_\_

Please fill out in block capitals!

| Pos. | Crowd | Decontaminated object | Serial number | Description / Comment |
|------|-------|-----------------------|---------------|-----------------------|
| 1    |       |                       |               |                       |
| 2    |       |                       |               |                       |
| 3    |       |                       |               |                       |
| 4    |       |                       |               |                       |

### Have the parts listed above been exposed to any of the following substances?

Health endangering watery solutions, buffers, acids, alkalis: ☐ Yes ☐ No

Potentially infectious agents: ☐ Yes ☐ No

Organic reagents and solvent: ☐ Yes ☐ No

Radioactive substances: ..... ☐ α.. ☐ β.. ☐ γ ..... ☐ Yes ☐ No

Health endangering proteins: ☐ Yes ☐ No

DNA: ☐ Yes ☐ No

These substances have reached the equipment/assembly? ☐ Yes ☐ No

Which one, if yes:

\_\_\_\_\_

Description of the measures for the decontamination of the listed parts:

\_\_\_\_\_  
\_\_\_\_\_

I confirm the proper decontamination:

Company/Dept. \_\_\_\_\_ Place and Date: \_\_\_\_\_

Signature of the authorized person n: \_\_\_\_\_

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P/N 30134185 D © 2017 OHAUS Corporation, all rights reserved / todos los derechos reservados / tous droits
Ohaus Corporation