
## Controls – Ranger 3000


| Button                                       | Primary Function<br>(Short Press <sup>1</sup> )                                                                                                       | Secondary Function<br>(Long Press <sup>2</sup> )                                          | Menu Function<br>(Short Press)                                                                                                                                                               |
|----------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div>On/Zero<br/>Off</div> <div>Yes</div>    | <ul style="list-style-type: none"> <li>• Turns scale on</li> <li>• Sets display to zero</li> </ul>                                                    | <ul style="list-style-type: none"> <li>• Turns the scale off</li> </ul>                   | <ul style="list-style-type: none"> <li>• Selects/accepts displayed setting</li> </ul>                                                                                                        |
| <div>Print<br/>Units</div> <div>No</div>     | <ul style="list-style-type: none"> <li>• Sends the current value to the selected COM ports if AUTOPRINT is set to Off.</li> </ul>                     | <ul style="list-style-type: none"> <li>• Toggles through active weighing units</li> </ul> | <ul style="list-style-type: none"> <li>• Advances to the next menu or menu item.</li> <li>• Rejects the current setting on the display and advances to the next available setting</li> </ul> |
| <div>Function<br/>Mode</div> <div>Back</div> | <ul style="list-style-type: none"> <li>• Initiates an application mode</li> </ul>                                                                     | <ul style="list-style-type: none"> <li>• Allows changing the application mode</li> </ul>  | <ul style="list-style-type: none"> <li>• Returns to previous menu item</li> </ul>                                                                                                            |
| <div>M+<br/>Menu</div> <div>Exit</div>       | <ul style="list-style-type: none"> <li>• Accumulates the weight and displays the accumulated information at 0 load.</li> </ul>                        | <ul style="list-style-type: none"> <li>• Enters User menu</li> </ul>                      | <ul style="list-style-type: none"> <li>• Quickly exit User Menu.</li> <li>• Aborts calibration in progress</li> </ul>                                                                        |
| <div>Tare</div>                              | <ul style="list-style-type: none"> <li>• Enter/clear a tare value. Clears the accumulation when the accumulation information is displayed.</li> </ul> |                                                                                           |                                                                                                                                                                                              |

**Notes:** <sup>1</sup> Short Press: Press less than 1 second.


<sup>2</sup> Long Press: Press and hold for more than 2 seconds.

## Controls – Ranger Count 3000

| Button                       | Primary Function<br>(Short Press)                                                                        | Secondary Function<br>(Long Press)                                                                              | Menu Function<br>(Short Press)                                                                                                 | Library Function<br>(Short Press)                                                                                                 |
|------------------------------|----------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>On/Zero Off</b><br>Yes    | <ul style="list-style-type: none"> <li>Turn scale on</li> <li>Zero the display</li> </ul>                | <ul style="list-style-type: none"> <li>Turn scale off.</li> </ul>                                               | <ul style="list-style-type: none"> <li>Accept the current menu or setting.</li> </ul>                                          | <ul style="list-style-type: none"> <li>Accept the current setting.</li> </ul>                                                     |
| <b>Print Units</b><br>No     | <ul style="list-style-type: none"> <li>Send the displayed value to the COM port</li> </ul>               | <ul style="list-style-type: none"> <li>Change the weighing unit.</li> </ul>                                     | <ul style="list-style-type: none"> <li>Advance to the next menu or setting.</li> <li>Increment the displayed value.</li> </ul> | <ul style="list-style-type: none"> <li>Advance to the next library or setting.</li> <li>Increment the displayed value.</li> </ul> |
| <b>Target</b><br>Back        | <ul style="list-style-type: none"> <li>Initiate the function of the current application mode.</li> </ul> | <ul style="list-style-type: none"> <li>Switch between Check Weigh, Check Count and Off.</li> </ul>              | <ul style="list-style-type: none"> <li>Previous menu</li> <li>Decrement the displayed value.</li> </ul>                        | <ul style="list-style-type: none"> <li>Previous library</li> <li>Decrement the displayed value.</li> </ul>                        |
| <b>M+ Menu</b><br>Exit       | <ul style="list-style-type: none"> <li>Accumulate the weight or pieces.</li> </ul>                       | <ul style="list-style-type: none"> <li>Enter user menu.</li> <li>View the Audit Trail event counters</li> </ul> | <ul style="list-style-type: none"> <li>Exit the menu.</li> <li>Abort the calibration</li> </ul>                                | <ul style="list-style-type: none"> <li>Exit the library.</li> </ul>                                                               |
| <b>ID</b>                    | <ul style="list-style-type: none"> <li>ID number input</li> </ul>                                        |                                                                                                                 |                                                                                                                                |                                                                                                                                   |
| <b>APW</b>                   | <ul style="list-style-type: none"> <li>Display/Store APW</li> </ul>                                      |                                                                                                                 |                                                                                                                                |                                                                                                                                   |
| 7 8 9<br>4 5 6<br>1 2 3<br>0 | <ul style="list-style-type: none"> <li>Enter numeric values on the display.</li> </ul>                   |                                                                                                                 |                                                                                                                                |                                                                                                                                   |
| <b>1 ΔΔ 2</b>                | <ul style="list-style-type: none"> <li>Enter dot (.) on the display.</li> </ul>                          | <ul style="list-style-type: none"> <li>Switch between scale 1 and scale 2</li> </ul>                            |                                                                                                                                |                                                                                                                                   |
| <b>C +/-</b>                 | <ul style="list-style-type: none"> <li>Clear last character from the display.</li> </ul>                 | <ul style="list-style-type: none"> <li>Switch between positive and negative value</li> </ul>                    |                                                                                                                                |                                                                                                                                   |
| <b>Tare</b>                  | <ul style="list-style-type: none"> <li>Perform a tare operation.</li> </ul>                              |                                                                                                                 |                                                                                                                                |                                                                                                                                   |
| <b>Sample</b>                | <ul style="list-style-type: none"> <li>Display/Start APW</li> </ul>                                      |                                                                                                                 |                                                                                                                                |                                                                                                                                   |

## Controls – Valor 7000

| Button                                                                              | Primary Function<br>(Short Press <sup>1</sup> )                                                                                                     | Secondary Function<br>(Long Press <sup>2</sup> )                                        | Menu Function<br>(Short Press)                                                       |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| <div>On/Zero</div> <div>Off</div> <div>Yes</div>                                    | <ul style="list-style-type: none"> <li>Turns scale on</li> <li>Sets display to zero</li> </ul>                                                      | <ul style="list-style-type: none"> <li>Turns the scale off</li> </ul>                   | <ul style="list-style-type: none"> <li>Selects/accepts displayed setting</li> </ul>  |
| <div>Print</div> <div>Units</div> <div>No</div>                                     | <ul style="list-style-type: none"> <li>Sends the current value to the selected COM ports if AUTOPRINT is set to Off.</li> </ul>                     | <ul style="list-style-type: none"> <li>Toggles through active weighing units</li> </ul> | <ul style="list-style-type: none"> <li>Toggles through available settings</li> </ul> |
| <div>Function</div> <div>Mode</div> <div>Back</div>                                 | <ul style="list-style-type: none"> <li>Selects function setting</li> </ul>                                                                          | <ul style="list-style-type: none"> <li>Selects active Mode</li> </ul>                   | <ul style="list-style-type: none"> <li>Returns to previous settings</li> </ul>       |
| <div>M+</div> <div>Menu</div> <div>Exit</div>                                       | <ul style="list-style-type: none"> <li>Accumulates the weight and displays the accumulated information at 0 load.</li> </ul>                        | <ul style="list-style-type: none"> <li>Enters User Menu</li> </ul>                      | <ul style="list-style-type: none"> <li>Quickly exit User Menu</li> </ul>             |
| <div>Tare</div>                                                                     | <ul style="list-style-type: none"> <li>Enter/clear a tare value. Clears the accumulation when the accumulation information is displayed.</li> </ul> |                                                                                         |                                                                                      |
|  | <ul style="list-style-type: none"> <li>The IR Sensor can be programmed to act as "touchless" button.</li> </ul>                                     |                                                                                         |                                                                                      |

**Notes:** <sup>1</sup> Short Press: Press less than 1 second.

<sup>2</sup> Long Press: Press and hold for more than 2 seconds.

Ohaus Corporation, 7 Campus Drive, Suite 310, Parsippany, New Jersey, 07054, USA [www.ohaus.com](http://www.ohaus.com)

**Declaration of conformity** We, Ohaus Corporation, declare under our sole responsibility that the balance models listed below marked with "CE" – are in conformity with the directives and standards mentioned.

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
For non-automatic weighing instruments used in an Article 1, 2.(a) application, additional metrological marking according to Annex IV of Council directive 2009/23/EC must be attached to the instrument.

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
Issued on: 2012-12-04

Ted Xia  
President  
Ohaus Corporation  
Parsippany, NJ USA

Robert Hansen  
Compliance Manager  
Ohaus Corporation  
Parsippany, NJ USA

\* 3 0 0 3 7 4 6 1 \*