


Reagent Bottle Narrow Mouth with Plastic Stopper Clear Glass



Part No	Size (ml)	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	F (mm)	G (mm)	H (mm)	I (mm)
LG1401-PS-60	60	44	61	81	89	36	22	14	16	21
LG1401-PS-125	125	57	77	100	103	36	22	14	16	21
LG1401-PS-250	250	69	102	130	140	40	25	18	20	25
LG1401-PS-500	500	80	127	155	160	48	31	22	23	29
LG1401-PS-1LTR	1000	110	150	184	192	52	31	22	25	32



Disclaimer: Dimensions given can vary slightly from manufacturer to manufacturer and are subject to change without notice. Dimensions should be used as a guide only. Wiltronics bears no responsibility for the accuracy of the dimensions given here.