



Reagent Bottle Narrow Mouth with Plastic Stopper Amber Glass


LG1402-PS-1LTRA
LG1402-PS-500A
LG1402-PS-250A
LG1402-PS-125A
LG1402-PS-60A

Part No	Size (ml)	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	F (mm)	G (mm)	H (mm)	I (mm)
LG1402- PS-60A	60	44	60	84	89	36	22	14	16	21
LG1402- PS-125A	125	54	77	102	109	36	22	14	16	21
LG1402- PS-250A	250	69	105	132	141	39	27	18	20	24
LG1402- PS-500A	500	80	127	153	160	49	31	22	23	29
LG1402- PS-1LTRA	1000	110	150	182	191	52	31	22	25	32



Scale: 2:1



Disclaimer: Dimensions given can vary slightly from manufacturer to manufacturer and are subject to change without notice. Dimensions should be used as a guide only. Wiltronics bears no responsibility for the accuracy of the dimensions given here.

www.wiltronics.com.au

Wiltronics Research Pty. Ltd. ABN 26 052 173 154
Unit 4 Cnr Ring Rd & Sturt St, Ballarat Victoria, 3350 | P.O Box 4043, Alfredton, 3350
sales@wiltronics.com.au | Phone: (03) 5334 2513 | Fax: (03) 5334 1845