


Reagent Bottle Narrow Mouth - Amber Glass



Part No	Size (ml)	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	F (mm)	G (mm)	H (mm)	I (mm)
LG1402-30A	30	39	46	64	74	30	27	20	14	17
LG1402-60A	60	46	59	81	92	36	30	22	17	21
LG1402-125A	125	57	70	90	111	37	31	23	17	22
LG1402-250A	250	69	95	123	140	40	36	27	20	22
LG1402-500A	500	80	127	157	165	47	36	27	23	29
LG1402-1LTRA	1000	105	150	182	196	50	37	29	24	31
LG1402-2.5LTRA	2500	143	255	275	292	69	56	45	38	45
LG1402-5LTRA	5000	170	275	320	340	83	72	60	44	57
LG1402-10LTRA	10000	210	330	400	420	83	72	60	44	57



Disclaimer: Dimensions given can vary slightly from manufacturer to manufacturer and are subject to change without notice. Dimensions should be used as a guide only. Wiltronics bears no responsibility for the accuracy of the dimensions given here.

www.wiltronics.com.au

Wiltronics Research Pty. Ltd. ABN 26 052 173 154
 Unit 4 Cnr Ring Rd & Sturt St, Ballarat Victoria, 3350 | P.O Box 4043, Alfredton, 3350
 sales@wiltronics.com.au | Phone: (03) 5334 2513 | Fax: (03) 5334 1845