


## Reagent Bottle Wide Mouth - Clear Glass



Part No	Size (ml)	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	F (mm)	G (mm)	H (mm)	I (mm)
LG1403-30	30	40	48	65	74	38	25	18	21	24
LG1403-60	60	48	62	80	90	43	27	19	24	28
LG1403-125	125	57	78	105	116	50	30	23	31	36
LG1403-250	250	69	98	125	138	58	36	27	38	42
LG1403-500	500	85	123	152	169	76	47	37	40	53
LG1403-1LTR	1000	110	158	190	203	78	42	31	49	55
LG1403-2.5LTR	2500	140	215	270	295	113	60	50	68	77



Scale: 2:1

**Disclaimer:** Dimensions given can vary slightly from manufacturer to manufacturer and are subject to change without notice. Dimensions should be used as a guide only. Wiltronics bears no responsibility for the accuracy of the dimensions given here.

[www.wiltronics.com.au](http://www.wiltronics.com.au)

Wiltronics Research Pty. Ltd. ABN 26 052 173 154  
 Unit 4 Cnr Ring Rd & Sturt St, Ballarat Victoria, 3350 | P.O Box 4043, Alfredton, 3350  
 sales@wiltronics.com.au | Phone: (03) 5334 2513 | Fax: (03) 5334 1845