

Reagent Bottle Blue GL45 Screw Cap - Clear Glass



LG1407-2000



LG1407-1000



LG1407-500




LG1407-250




LG1407-100

Part No	Size (ml)	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)
LG1407-100	100	57	66	103	51	27
LG1407-250	250	70	103	144	51	27
LG1407-500	500	88	140	178	51	27
LG1407-1000	1000	99	195	229	51	27
LG1407-2000	2000	137	230	262	51	27



Scale: 2:1
GL45 standard cap.



Disclaimer: Dimensions given can vary slightly from manufacturer to manufacturer and are subject to change without notice. Dimensions should be used as a guide only. Wiltronics bears no responsibility for the accuracy of the dimensions given here.