


Globe Type Separatory Funnels with Ground Glass Stopper and Stopcock

- Useful for separating liquids
- Ground glass joint at top
- Ground-in glass stopcock and stopper

These globe type separatory funnels feature a ground glass mouth and ground-in stopper and stopcock. They can be used to separate two immiscible liquids of different densities. The denser of the two liquids will move to the bottom of the globe of the separatory funnel and may then be removed via the stopcock.



Dimensions

Part No.	Size (ml)	A (mm)	B (mm)	C (mm)	D (mm)
LG1542-125	125	75	300	155	20
LG1542-250	250	85	340	165	20
LG1542-500	500	105	400	178	24

