

Glass Dropping Funnels with PTFE Valve Stopcock

- PTFE needle valve stopcock for controlling liquid flow
- 19/26 ground-in glass joint

These glass cylindrical dropping funnels are ideal for transferring fluids in the lab. They feature PTFE needle valve stopcocks to control the flow of liquid. They also have a ground-in glass sockets and cones.



Dimensions

Part No.	Capacity (ml)	Socket	Cone
LG5101-G-50C	50	19/26	19/26
LG5101-G-100C	100	19/26	19/26

