

610g Dial-O-Gram Balance

Features

- Dial operation for speed and convenience
- Tare capacity for fast repetitive weighings
- Notched and tiered beams with center reading
- A removable Stainless Steel platform
- Spring loaded zero adjust compensator



ME4163 - Ohaus Model No. 1610-00

Capacity (metric)	610 g / 2610 g (Capacity of 2610g requires purchase of optional 707-00 attachment weight set)
Readability (metric)	0.1 g / 0.1 g
Additional tare capacity (g)	225
Dial calibration (g)	10 x 0.1
Front beam calibration (g)	500 x 100
Rear beam calibration (g)	100 x 10
Platform Construction	Stainless Steel pan (removable)
Platform Size (Dia x Depth) (in/cm)	6 x 0.8 / 15.2 x 1.9
Net Wt. (lb/kg)	7 / 3.2
Shipping dimensions (WxDxH) (in/cm)	22 x 10 x 8 / 55.9 x 25.4 x 20.4
Shipping Wt. (lb/kg):	8 / 3.6

Distributed by

Wiltronics Research Pty. Ltd.

Unit 4 Cnr Ring Rd & Sturt Street Ballarat, Victoria 3350

Phone: (03) 5334 2513 Fax: (03) 5334 1845

Email: sales@wiltronics.com.au Web: www.wiltronics.com.au