

Ohaus Dissolved Oxygen Measurement Kit 0-19.9mg/L ME4355

Accurate Electrochemistry Measurement at Your Fingertips

- **Starter pen meters provide what should be expected of a small, economical pen meter: simple, fast, straightforward, no-frills operation time and time again**
- **With durable ABS housing, a protective sensor cap, and automatic shutdown feature that preserves battery life, starter pen meters can endure consistent use in rugged environments**
- **Equipped with a wrist strap that prevents unintentional dropping, starter pen meters feature an IP67 waterproof design that can prevent water damage when dropped into a liquid**



Specifications

Application	Dissolved oxygen (DO), temperature measurement
Display	Liquid crystal display (LCD) with secondary information line (on select models)
Operation	4 x 1.5V batteries (included)
Communication	N/A
Construction	Polargraphic electrode, uses membranes
Design Features	Easily replaceable electrodes, automatic shut off feature
Measurement Range	0.0-19.9mg/L
Calibration	100%, 0%
Item No.	30073986
Supplied With	Membrane caps (3), electrolyte solution and wrist strap
Accessories	Replacement electrode (item no. 30222083) Membrane kit (item no. 30222084)
Price	\$565.00