

Porcelain Low Form Crucibles with Lids

- Useful for heating applications in the lab
- Includes lid

Low form (squat form) porcelain crucibles with lids. Ideal labware for heating in science lab experiments. Pair these crucibles with a pipe clay triangle, tripod and Bunsen burner for a great heating setup.



Dimensions

Product No.	Size (ml)	A (mm)	B (mm)	C (mm)
LC1260-5	5	26	21	17
LC1260-15	15	22	22	42
LC1260-18	18	47	25	22
LC1260-25	25	49	27	25
LC1260-40	40	54	34	26
LC1260-50	50	59	37	30

