


Porcelain Tall Form Crucibles with Lids

- Useful for heating applications in the lab
- Includes lid

Tall form porcelain crucibles with lids. Ideal labware for heating in science lab experiments. Pair these crucibles with a pipe clay triangle, tripod and Bunsen burner for a great heating setup.



Dimensions

Product No.	Size (ml)	A (mm)	B (mm)	C (mm)
LC1262-20	20	38	41	24
LC1262-25	25	38	41	25
LC1262-30	30	39	45	25
LC1262-35	35	41	49	26
LC1262-40	40	42	52	26
LC1262-50	50	43	56	28

