

Porcelain Buchner Funnel Filters

- Made from acid- and alkali-resistant porcelain
- Fully autoclavable
- Designed to be used with a filtering flask and vacuum to speed up filtering process

These porcelain Buchner funnels are designed for filtration purposes in the lab. They are acid & alkali resistant. Buchner funnels are used in combination with a filtering flask so that vacuum suction can be applied, speeding up the filtration process considerably.



Dimensions

Product No.	Paper Dia. (mm)	A (mm)	B (mm)	C (mm)	D (mm)
LC1510-55	50	55	105	50	30
LC1510-70	70	75	115	55	33
LC1510-80	70/80	80	130	65	40
LC1510-90	90	90	140	70	40
LC1510-110	110	110	170	85	50
LC1510-125	125	125	175	80	60
LC1510-150	150	150	210	95	60

