

Glass Mortar & Pestles

- Spout for easy pouring
- Wall thickness of over 5mm for sturdiness

Fantastic for laboratory use, these heavy-duty glass mortar & pestles feature a wall thickness of over 5mm. They are ideal for grinding solids into smaller pieces or a fine powder. A great tool to have for the lab.

Dimensions

Product No.	Size (ml)	A (mm)	B (mm)	C (mm)	D (mm)
LG1383-60	60	60	41	43	105
LG1383-75	75	75	49	60	105
LG1383-90	90	90	59	70	117
LG1383-120	120	120	75	88	120
LG1383-150	150	150	82	104	170

