

## Straight Thistle Funnels

- Design allows for safe handling of small volumes of liquid
- Available with straight stem or in 1 or 2 bulb variants
- Made from high quality glass

These thistle funnels are designed for safe liquid handling. They may be passed through a bung or clamped in place. Thistle funnels are useful when a small volume of liquid must be added to another piece of labware and other funnels have stems that are not long enough.



### Dimensions

Part No.	Type	Length (mm)	Top Diameter (mm)	Stem Diameter (mm)
LG1521-1	Straight	350	40	7-8
LG1523-1	1 Bulb	350	40	7-8
LG1524-1	2 Bulb	350	40	7-8

