

Reagent Storage Square Jar - Polypropylene

LP1405PP-500, LP1405PP-1000



Description

Polypropylene square reagent (media) storage jar used to store liquid or powder chemicals for science experiments in the lab. Storage space is valuable in any laboratory. These storage jars are essential for any science lab or classroom requiring an inexpensive and robust way to safely keep chemicals.

Part No.	Size (ml)	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)
LP1405PP-500	500	80	95	118	60	16
LP1405PP-1000	1000	105	115	140	82	17

