

## Mini Vortex Mixer (Fixed Speed)

**OHAUS VXMNFS**

**LQ1510FS**

- **Sturdy design**
- **Fixed high speed mixing**
- **Touch mode operation**

### Description

Built sturdy to provide stable and reliable vortexing action. Starts mixing when the cup head is pressed down. Speed is fixed at full RPM to provide vigorous vortexing of samples.

### Operating Conditions

Unit can be run in conditions from 4 to 40°C, 20% to 85% relative humidity, non-condensing.

### Includes

Cup Head, 7.6 cm Head, 7.6 cm Rubber Head Cover

### Specifications

**Speed Range 230V (RPM)** 2500

**Orbit (mm)** 4.9

**Controls** None

**Duty Rating** Intermittent duty

**Dimensions (L x W x H cm)** 21.1 x 12.2 x 16.5

**Ship Weight (kg)** 4.5



### Dimensions

