

Magnetic Stirrer Ohaus Mini LQ2072

- Easy to clean, chemical resistant ceramic top plate
- Adjustable speed control
- Compact size



Description

The OHAUS Mini Magnetic Stirrer is a rugged, compact unit that can stir up to 1000 ml of liquid. These stirrers are ideal for educational and other lab environments that require consistent stirring. It features a 102 x 102 mm, easy-to-clean ceramic top plate and valox polymer housing.

Specifications

Function	Lab Magnetic Stirring
Capacity	1000 mL
Top Plate Construction	Ceramic
Top Plate Size	102 x 102 mm
Control	Analog
Power	230V, 0.25A, 50/60Hz
Power Consumption	20 W
Safety Certification	CE; TUV
Speed Range	100 rpm to 1200 rpm
Working Environment	5°C to 50°C, 80%RH, non-condensing
Dimensions (L x W x H)	20.1 x 15.2 x 12.7 cm

