

Hotplate Stirrer Mini Ohaus LQ2074

- Adjustable speed magnetic stirrer
- Adjustable temperature hotplate
- Easy to clean, chemical resistant ceramic hotplate
- Cool touch valox polymer housing
- Compact size



Description

The OHAUS Mini Hotplate with Magnetic Stirrer is a rugged, compact units that can heat and stir up to 1000 ml of liquid. These hotplate/stirrers are ideal for educational and other lab environments that require heating and consistent stirring. It features a 102 x 102 mm, easy-to-clean ceramic top plate and cool-touch housing.

Specifications

Function	Heating and Magnetic Stirring
Capacity	1000 mL
Hot Plate Construction	Ceramic
Hot Plate Size	102 x 102 mm
Temperature Range	Ambient temperature to 500°C
Control	Analog Adjustable Temperature and Adjustable Stirring Speed
Power	230V, 0.25A, 50/60Hz
Power Consumption	200 W
Safety Certification	CE; TUV
Speed Range	100 rpm to 1200 rpm
Working Environment	5°C to 50°C, 80%RH, non-condensing
Dimensions (L x W x H)	20.1 x 15.2 x 12.7 cm

