

## **Technical Data Sheet**

## Vitals 1209 Vinyl Powder Free Glove Item# VITUP1209X

SAP Material Code & Size	SAP Material Code	Product Reference Code	Size
	VITUP1209B-0100100	VITUP1209B	Small (S)
	VITUP1209C-0100100	VITUP1209C	Medium (M)
	VITUP1209D-0100100	VITUP1209D	Large (L)
	VITUP1209E-0100100	VITUP1209E	Extra Large (XL)
Country of Origin	China		
Intended Use	Used to avoid cross-contamination. Applicable to single-use gloves that are worn once and then discarded.		
Length & Width	Size	Total Length (mm)	Palm Width (mm)
	Small (S)	Total Longth (Inny	85
	Medium (M)		95
	Large (L)	Min 235	105
	Extra Large (XL)		115
Thickness	Min 0.05 mm (Finger) Min 0.06 mm (Palm)		
Material	PVC		
Color	Clear		
Texture	Smooth		
Elongation	Before Aging $\ge$ 300% After Aging $\ge$ 300%		
Tensile Strength	Before Aging $\ge$ 11 MPa After Aging $\ge$ 11 MPa		
Powder Residue	≤ 2.0 mg per glove		
Sterility	Non-sterile		
Usage	Single-use only, disposable		
Packaging	100 gloves per box,10 boxes per case / 1000 gloves per case		
Shelf Life	3 years		
Recommended Storage Conditions	Store in a cool and dry place. Avoid excessive heat (40°C, 104°F). Shield open box from direct sun and fluorescent lighting.		

This technical data sheet is issued by A.R. Medicom Inc. (Asia) Ltd.