

**TECHNICAL DATA SHEET**  
**Nitrile, Non-Sterile Medical Examination Glove**  
**which meet ASTM D 6319-19**

**Intended Use**

Non-sterile, single use, medical examination gloves are intended to prevent transmission of a wide variety of diseases to both patients and health care personnel. Specifically, gloves are intended to act as a barrier between the health care provider and patients to prevent contamination between patient and caregiver.

Characteristic		Acceptable Limit	
Dimension & Tolerance <sup>(1)</sup>	Width (mm)	Extra-Small (XS)	60 - 80
		Small (S)	70 - 90
		Unisize (S-M)	75 - 95
		Medium (M)	85 - 105
		Large (L)	100 - 120
		X-Large (XL)	110 - 130
		XX-Large (XXL)	120 - 140
	Length (mm)	XS & S: ≥ 220 Others: ≥ 230	
Thickness (mm)	Palm	≥ 0,05	
	Finger		
Physical Requirements	Before ageing	Tensile Strength	≥ 14 MPa
		Ultimate Elongation	≥ 500%
	After ageing	Tensile Strength	≥ 14 MPa
		Ultimate Elongation	≥ 400%
Powder Residue	Powder Free Gloves	≤ 2,0 mg/glove	
	Powdered Gloves	≤ 10,0mg/dm <sup>2</sup>	
Description	Ambidextrous. Refer to Appendix 1 for color and packaging information.		

<sup>(1)</sup> Dimensions and Tolerances vary with manufacturers.

**Additional Information**

Storage Condition	Country of Origin
Store under cool, dry place. Avoid excessive heat (40°C, 104°F). Shield open box from direct sun and fluorescent light.	Malaysia or China or Thailand

**Appendix 1 List of Item Codes and Attributes**

## Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
1137	Medicom SafeTouch "Advanced Stretchy" Powder-Free Nitrile Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Country of origin: Malaysia	Non-Sterile
1174	Medicom SafeTouch "Advanced Platinum" Powder-Free Nitrile Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	White	Country of origin: Malaysia	Non-Sterile
1175	Medicom SafeTouch "Advanced Slim" Powder-Free Nitrile Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Country of origin: Malaysia	Non-Sterile
1176	Medicom SafeTouch "Advanced Rejuvenate Powder-Free Nitrile Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	White	Country of origin: Malaysia	Non-Sterile

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
1177	Medicom SafeTouch "Advanced Form" Powder-Free Nitrile Medical Examination Gloves	200 Gloves/Box, 10 Boxes/Case, 2000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Country of origin: Malaysia	Non-Sterile
1185	Medicom SafeTouch "True Fit Thin" Powder-Free Nitrile Medical Examination Gloves	300 Gloves/Box, 8 Boxes/Case, 2400 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Country of origin: Malaysia	Non-Sterile
1185EX	Medicom Safe Basics "True Fit Thin" Powder-Free Nitrile Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Country of origin: Malaysia	Non-Sterile
1186	Medicom SafeTouch "True Fit Thin" Powder-Free Nitrile Medical Examination Gloves	300 Gloves/Box, 8 Boxes/Case, 2400 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large	Pink	Country of origin: Malaysia	Non-Sterile
1196	Medicom SafeTouch "Advanced Comfort" Powder-Free Nitrile Medical Examination Gloves	300 Gloves/Box, 10 Boxes/Case, 3000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	White	Country of origin: China	Non-Sterile
1199	Medicom SafeTouch "Advanced Rejuvenate Powder-Free Nitrile Medical Examination Gloves	200 Gloves/Box, 10 Boxes/Case, 2000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large	White	Country of origin: Malaysia	Non-Sterile

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
1217	Medicom SafeTouch "Advanced Extend+" Powder-Free Nitrile Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large	Pink	Country of origin: Malaysia	Non-Sterile
1220	Medicom SafeTouch "Advanced™ Force" Powder-Free Nitrile Medical Examination Gloves	100 Gloves/Box, 10 Boxes/Case, 1000 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Blue	Country of origin: Malaysia	Non-Sterile
1227	Medicom SafeTouch "SilkSkin" Powder-Free Nitrile Medical Examination Gloves	150 Gloves/Box, 10 Boxes/Case, 1500 Gloves/Case Gloves are in a cardboard. Case made with suitable corrugated fiberboard to withstand transportation and storage	A = X-small, B = Small, C = Medium, D = Large, E = X-Large	Indigo	Country of origin: Malaysia	Non-Sterile