

Newton Meter NaRiKa GN-1

ME4150



**Measuring force became
so quick and precise!**

Digital display meter instead of traditional balance scales

Features

- Quick and precise measuring of force
- Digital display
- Zero point adjustment function
- Hold function
- Measures both pushing and pulling forces
- Two modes: Newton (N) / Weight (g)

Specifications

- Measuring Units: Newton (N) / Weight (g)
- Measurement Range: $\pm 19.99 \text{ N}$ / $\pm 1,999 \text{ g}$
- Resolution: 0.01 N / 1 g