

Data Sheet

MS1010

Stereo Binocular Microscope x20 LED Illuminations

- Perfect for primary and junior secondary educational use
- LED Light Sources
- Portable – battery powered
- Great for outdoor use

The MS1010 Stereo Microscope is a small light weight microscope, ideal for primary and junior secondary use. It is battery operated (2x AA batteries) making it extremely safe and portable for use in the classroom and out in the field. Great for viewing plant materials (leaves, flowers etc.), feathers, insects, coins and stamps.

Specifications

Magnification	20X
Viewing Head	Upright Binocular Viewing Head
Eyepiece	Wide Field (WF) 10X
Objective	Objective 2X
Inter-pupillary distance	55 - 75mm
Working Distance	70mm
Illumination	Built-in Incident LED Light
Power	2X AA Batteries (not included)
Weight	600g
Overall Height	230mm



Page 1 of 1



**Your Science & Technology
Supply Specialists**

www.wiltronics.com.au



Wiltronics Research Pty. Ltd. ABN 26 052 173 154
5-7 Ring Road, Alfredton Victoria 3350 | P.O Box 4043, Alfredton, 3350
sales@wiltronics.com.au | Phone: (03) 5334 2513 | Fax: (03) 5334 1845