



# Monocular Biological Microscope

40x - 400x, LED Illumination

**MS1020A**

- **Magnification 40x - 400x**
- **Coarse & Fine Focus Adjustment**
- **LED Transmissive Light Source**
- **LED Lamp Adjustable Brightness**

THE MS1020A is a very popular entry level biological microscope used for observing biological slices, some biological cells, transparent or translucent objects and powders, fine particles and other objects



## Product Specifications

<b>Head</b>	Monocular head, 30° inclined, 360° rotatable
<b>Eyepiece</b>	WF 10x with pointer
<b>Magnification</b>	40x-1000x
<b>Objective</b>	Achromatic 4x Achromatic 10x Achromatic 40x, spring
<b>Nosepiece</b>	Triple revolving nosepiece
<b>Stage</b>	Plain stage with clips, size 100mm x 90mm
<b>Condenser</b>	Single lens condenser with Disc diaphragm
<b>Focusing</b>	Coaxial Coarse & Fine Focusing Knobs
<b>Lighting</b>	LED Lamp adjustable output
<b>Power Supply</b>	240V IEC Lead

