



Microscope Binocular Monocular 40x - 400x, LED Illumination

- **Magnification 40x - 400x**
- **Adjustable brightness**
- **LED transmissive light source**
- **Coarse & Fine Focus adjustment**

THE MS1020A is a very popular entry level biological microscope used for observing biological slices, some biological cells, transparent or translucent objects and powders, fine particles and other objects

Applications

Ideal for educational use in biology for viewing prepared slides of biological slices, biological cells, bacteria and living tissue cultures, fine particles and other objects



MS1020A

Product Specifications

Head	Monocular Head, 30 degrees inclined, 360 Degree Rotatable
Eyepiece	WF10X with pointer
Nosepiece	Triple Revolving Nosepiece
Objective	Achromatic 4x Achromatic Objectives 10x Achromatic Objectives 40x, spring
Stage	Plain Stage with Clips, Size 100mm x 90mm
Condenser	Single lens condenser with Disc diaphragm
Focus	Coaxial Coarse & Fine Focusing Knobs
Light	LED Lamp, adjustable output
Power Supply	240VAC to 5VDC Plug Pack Adaptor

