



# Monocular Biological Microscope

40x - 1000x, LED Illumination

**MS1035A**

- **Magnification 40x - 1000x**
- **Wide Field Plan WF10X**
- **Quadruple Nosepiece**
- **LED Lamp Adjustable Brightness**

The Monocular Biological Microscope, 40x - 1000x, LED Illumination, MS1035A is a great entry-level biological microscope, especially with a limited budget. It features an integrated X-Y mechanical stage and an Abbe N.A. 1.25 condenser.



## Product Specifications

<b>Head</b>	Monocular Head, 30° Inclined, 360° Rotatable
<b>Eyepiece</b>	Wide Field Plan WF10X
<b>Magnification</b>	40x-1000x
<b>Objective</b>	Achromatic 4x Achromatic 10x Achromatic 40x (s) Achromatic 100x (S, Oil)
<b>Nosepiece</b>	Quadruple Nosepiece
<b>Stage</b>	Double Layers Mechanical Stage 140 x 155 mm, moving range 75 x 50mm
<b>Condenser</b>	Abbe Condenser N.A. 1.25, $\Phi 2\text{-}\Phi 30\text{mm}$ Iris Diaphragm, with $\Phi 32$ filter holder
<b>Focusing</b>	Coaxial Coarse & Fine Focusing Adjustable Mechanism 25 mm, Precision 0.002 mm, with Tension Adjustment, with Safety Stop Ring
<b>Lighting</b>	LED Lamp 1W, brightness adjustable
<b>Power Supply</b>	240VAC to 5V 1A Power Plug Pack

