



MS1044A-DM WIFI Digital Microscope



**Built-in WIFI Camera
2MP. With Pad Holder**

**Plan Eyepiece
10x/18mm With Eye Guard**

**0.2W LED(5700-6500K),
Intensity Adjustable**



ElectroTech
Instruments

Digital WIFI Microscope

Binocular Head, 360° Rotatable

Built-in WIFI Camera
Quadruple Nosepiece
Adjustable LED
Adapted Power Supply
Wi-Fi Display Connection

*Screen Pictured is not Included
** Pictured with a Wi-Fi Connected
Android Tablet



MS1044A-DM

Product Specifications

Head	Binocular Head 30° Inclined, 360° Rotatable, Inter-pupillary Distance 50-75mm, Built-in WIFI Camera, 2MP. With Pad Holder
Eyepiece	Plan Eyepiece 10x/18mm With Eye Guard
Nosepiece	Built-in Quadruple Nosepiece
Objective	Achromatic Objectives 4x/0.10(Black) Achromatic Objectives 10x/0.25(Black) Achromatic Objectives 40xS/0.65(Black) Achromatic Objectives 100xS/1.25(Oil)(Black)
Stage	115x125mm Double Layers Mechanical Stage, With N.A.1.25 Abbe Condenser, With Iris Diaphragm & Filter Holder, Moving Range 22mm-71mm
Focus	Coaxial Coarse & Fine Focus, With Coarse Adjustment Stop Screw And Tightness Adjustment, Focusing Range 13mm With Clip Stage, (7mm With Mechanical Stage), Fine Precision 0.004mm
Light	0.2W LED(5700-6500K), Intensity Adjustable
Power	12V DC Plug pack