



# Digital USB Microscope

50x - 1000x, LED Illumination

- **Magnification 50x - 1000x**
- **Universal Rotating Base**
- **USB Connection**
- **8 Built In LEDs**



**MS1312**

## Product Specifications

<b>Magnification</b>	50x-1000x
<b>Image Sensor</b>	1/6.5" Color CMOS Image Sensor
<b>Resolution</b>	640x480, 30Fps
<b>White Balance</b>	Auto
<b>Connection</b>	USB 2.0
<b>Illumination</b>	8 LEDs, Brightness Adjustable
<b>USB Support</b>	Abbe Condenser N.A. 1.25, $\Phi$ 2- $\Phi$ 30mm Iris Diaphragm, with $\Phi$ 32 filter holder
<b>Output</b>	BMP, JPEG
<b>SNR</b>	Max 42dB

