



MS1316

Digital Pen Microscope Stand with Adjustable Stage

- Easy setup, perfect for classroom
- Battery or USB powered
- USB cable included
- 11 LED transillumination with precise brightness adjustment
- Stage is adjustable on X and Y axes
- Design makes it suitable for holding a number of different pen microscopes

The Digital Pen Microscope Stand features a diversified holder design that makes the stand suitable for a range of different pen microscopes. The position of the stage is adjustable on the X & Y axis for optimum precision. The base of the stand boasts 11 LED transillumination with precise brightness adjustment. The stand is very simple to setup, making it perfect for a classroom environment. It comes with a USB cable which can power the stand, or alternatively it can be battery-powered.

