



Digital Pen Microscope

12MP, USB, on Stand

- **10x - 220x Magnification**
- **Direct HDMI Connection**
- **1080p Video Recording**
- **6 x LED Lights**

A microscope for use by students, teachers, hobbyists, collectors etc. Provides USB PC interface plus a direct output to a 1080p TV screen via HDMI. Super sharp optics offer up to 220x magnification with adjustable focus between 10 and 50mm. The 2.4" LCD screen allows it to be used standalone without a PC if required. Takes still shots up to 12mp in size for close up inspection. It even records 1080p HD video to internal SD card (or via USB PC connection). An inbuilt battery adds portability and caters for several hours use.



MS1325

Product Specifications

Magnification	10x to 220x
Display	2.4 Inch LCD
Photo Size	3M, 5M, 8M, 12M
Video Size	1080p, 720p, VGA
Video Output	HDMI, TV, USB
Frame Rate	Max 30 f/s under 600 Lux Brightness
Image Sensor	3 MP
Light Source	6 LED (adjustable by control wheel)
Power Source	Built-in Li-ion Battery or 5V DC-in
Storage	Micro-SD Card up to 32GB

