



Wi-Fi Digital Microscope

50x - 1000x, 2MP, USB

- **Optical Zoom 50x - 1000x**
- **HD Colour CMOS sensor**
- **Adjustable LED Brightness**
- **Quick Capture Button**

This Wi-Fi Digital Microscope has a 2 MP camera and can share a real time image to mobile phones and tablets. It has a 10 metre Wi-Fi range. This microscope also features adjustable LED brightness and 50x-1000x magnification. An easy, compact and convenient microscope with 3 hours of battery life.



MS1330

Product Specifications

| | |
|-----------------------|------------------------------------|
| Magnification | 50x-1000x |
| Lens Diameter | 8.0mm |
| Resolution | 640x480, 1920x1080 |
| Camera Pixel | 2.0MP |
| Battery | 3 hours Battery Life |
| WIFI Distance | WIFI 10m Range |
| Illumination | 8 LEDs |
| WIFI Support | Android, IOS, USB Support Win, Mac |
| Output | MJPEG/AVI |
| Color | Black |
| Focusing Range | 3~60mm |

